

MyWay+ Measurement Guide




Child name: _____

Date: _____

Location: _____

Therapist: _____

Notes: _____



- A** Thigh length (Back of knee to back of bottom)
- B** Hip Width (From either side of widest part of hips)
- C** Back height (From bottom to top of shoulder)
- D** Chest width (from either side of chest, nipple level)
- E** Lower leg length (From back of knee to underside of foot)

Frame Size	1	2	3
Age (approx.)	1-5	4-11	9-16
User Shoulder Height (to floor) (inches)	23 - 34"	34 - 47"	42 - 59"
User Saddle Height (inside leg) (inches)	9.5 - 16"	16 - 24"	21.5 - 32"
Max. User Weight (lb)	55 lbs.	110 lbs.	176 lbs.
Frame Weight (lb)	26.5 lbs	14.5 lbs.	21 lbs.
Upper Support Unit Weight (lb)	Size 1 frame & upper support unit do not separate	23 lbs.	24 lbs.
Frame Width (inches)	22.5"	26.5"	28.5"
Frame Length (inches)	33"	36.5"	46.5"

Harness Size	1	2	3	4	5	6
Age (approx.)	1-4	2-6	5-9	8-12	12-15	Up to 16
Waist at Navel Level (inches)	18 - 22"	19.5 - 23.5"	21.5 - 25.5"	23.5 - 27.5"	27.5 - 31.5"	31.5 - 35.5"
Max. User Weight (lb)	110 lbs.	110 lbs.	110 lbs.	176 lbs.	176 lbs.	176 lbs.
Color	Purple	Green	Blue	Red	Orange	Grey