

Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC		Phone	
Mark For				Email / Fax			

✓ OPTION	DESCRIPTION	TOTAL STANDARD PACKAGE	RETAIL
Base			\$

Leckey only compatible with casters, frame size, Leckey chair size, upholstery color, and accessory selections

☐ 801BS65	MC4000Ped Must select 801CW123 or 801CW124 Casters	\$ 3,684.00	\$ 2,525.00
☐ 801BS66	MC4000TiltPed Must select 801TT10 Add Tilt, and 801CW123 or 801CW124 Casters	\$ 5,133.00	\$ 2,525.00
☐ 801BS68	MC4020Ped Must select 801WQ12 or 801WQ15 20" Rear Wheels w/ 5" or 6" Casters	\$ 4,623.00	\$ 2,525.00
☐ 801BS67	MC4000Leckey Must select 801CW123 or 801CW124 Casters and Leckey Seat on page pg.3	\$ 4,133.00	\$ 1,322.00

Base Options	\$	\$
--------------	----	----

Tilt

Only add tilt if ordering the MC4000TiltPed

☐ 801TT10	Add Tilt	\$ 1,449.00
----------------------------------	----------	-------------

Base Type

Casters	\$	\$
---------	----	----

Only select casters if ordering a rolling base

☐ 801CW123	5" Casters (all four legs)	\$ 1,159.00	☐ 801CW124	6" Casters (all four legs)	\$ 1,159.00
-----------------------------------	----------------------------	-------------	-----------------------------------	----------------------------	-------------

Wheels	\$	\$
--------	----	----

Only select wheels if ordering the MC4020Ped

☐ 801WQ12	20" Rear Wheels with 5" Casters	\$ 2,098.00	☐ 801WQ15	20" Rear Wheels with 6" Casters	\$ 2,098.00
----------------------------------	---------------------------------	-------------	----------------------------------	---------------------------------	-------------

Frame Size	\$
------------	----

☐ 801W70	17" W x 17" D Frame (not compatible with MC4000Leckey with size 3 or 4 Leckey chair)	STD
☐ 801W71	17" W x 20" D Frame (compatible with MC4000Leckey with size 3 or 4 Leckey chair only)	\$ 414.00

Seat Depth	\$	\$
------------	----	----

*16.20 & 17.70" seat depth not compatible with tall & narrow back or Laterals. *16.90 & 17.70" seat depth requires seat depth adapters.

☐ 801SD109	10.70" Seat Depth	NC	☐ 801SD113	15.40" Seat Depth	NC
☐ 801SD110	11.80" Seat Depth	NC	☐ 801SD114	16.20" Seat Depth*	NC
☐ 801SD111	13.20" Seat Depth	NC	☐ 801SD115	16.90" Seat Depth*	NC
☐ 801SD112	14.30" Seat Depth	NC	☐ 801SD116	17.70" Seat Depth**	NC

Seat Cushion	\$	\$
--------------	----	----

☐ 5.5" Wide Commode Opening with 4" Slot	☐ 801CU140	No Gel	STD	☐ 801CU143	Gel	\$ 799.00
☐ 5.5" Wide Commode Opening with No Slot	☐ 801CU141	No Gel	\$ 145.00	☐ 801CU144	Gel	\$ 799.00
☐ No Commode Opening (solid seat)	☐ 801CU142	No Gel	\$ 164.00	☐ 801CU145	Gel	\$ 799.00
☐ Custom Cushion (select custom seat cushion in next section)						BELOW

Quick Release Pins

Quick release pins for two seat cushions requires one of the cushions to be no commode opening (solid seat).

☐ 801CU133	Add Quick Release Pins for 1 Seat Cushion (easily take seat cushion off frame for cleaning)	\$ 81.00
☐ 801CU146	Add Quick Release Pins for 2 Seat Cushions (easily switch between two seat cushions) (select two seat cushions)	\$ 163.00

✓ OPTION	DESCRIPTION							RETAIL	
Custom Seat Cushion					\$		\$		
Select style and size/gel option									
☐ Tear-Drop Open (TDOP)		1-Part	4" Opening x 3.0" Slot	☐ 801QC63	No Gel	\$ 642.00	☐ 801QC87	Gel	\$ 799.00
			5" Opening x 3.0" Slot	☐ 801QC64	No Gel	\$ 642.00	☐ 801QC88	Gel	\$ 799.00
			6" Opening x 3.5" Slot	☐ 801QC65	No Gel	\$ 642.00	☐ 801QC89	Gel	\$ 799.00
☐ Racetrack Open (RTOP)		1-Part	4" Opening x 3.0" Slot	☐ 801QC66	No Gel	\$ 642.00	☐ 801QC90	Gel	\$ 799.00
			5" Opening x 3.0" Slot	☐ 801QC67	No Gel	\$ 642.00	☐ 801QC91	Gel	\$ 799.00
			6" Opening x 3.5" Slot	☐ 801QC68	No Gel	\$ 642.00	☐ 801QC92	Gel	\$ 799.00
☐ Tear-Drop Closed (TDCP)		1-Part	4" Opening x No Slot	☐ 801QC69	No Gel	\$ 642.00	☐ 801QC93	Gel	\$ 799.00
			5" Opening x No Slot	☐ 801QC70	No Gel	\$ 642.00	☐ 801QC94	Gel	\$ 799.00
			6" Opening x No Slot	☐ 801QC71	No Gel	\$ 642.00	☐ 801QC95	Gel	\$ 799.00
☐ Racetrack Closed (RTCP)		1-Part	4" Opening x No Slot	☐ 801QC72	No Gel	\$ 642.00	☐ 801QC96	Gel	\$ 799.00
			5" Opening x No Slot	☐ 801QC73	No Gel	\$ 642.00	☐ 801QC97	Gel	\$ 799.00
			6" Opening x No Slot	☐ 801QC74	No Gel	\$ 642.00	☐ 801QC98	Gel	\$ 799.00
☐ Tear-Drop Side Cutout Left Side (TDSLPL)		1-Part	4" Opening x 3.0" Slot	☐ 801QC77	No Gel	\$ 642.00	☐ 801QC101	Gel	\$ 799.00
			5" Opening x 3.0" Slot	☐ 801QC75	No Gel	\$ 642.00	☐ 801QC99	Gel	\$ 799.00
			6" Opening x 3.5" Slot	☐ 801QC76	No Gel	\$ 642.00	☐ 801QC100	Gel	\$ 799.00
☐ Tear-Drop Side Cutout Right Side (TDSRP)		1-Part	4" Opening x 3.0" Slot	☐ 801QC80	No Gel	\$ 642.00	☐ 801QC104	Gel	\$ 799.00
			5" Opening x 3.0" Slot	☐ 801QC78	No Gel	\$ 642.00	☐ 801QC102	Gel	\$ 799.00
			6" Opening x 3.5" Slot	☐ 801QC79	No Gel	\$ 642.00	☐ 801QC103	Gel	\$ 799.00
☐ Straight U-Slot (SUSP)		1-Part	4" Opening x 4" Slot	☐ 801QC81	No Gel	\$ 642.00	☐ 801QC105	Gel	\$ 799.00
			5" Opening x 5" Slot	☐ 801QC82	No Gel	\$ 642.00	☐ 801QC106	Gel	\$ 799.00
			6" Opening x 6" Slot	☐ 801QC83	No Gel	\$ 642.00	☐ 801QC107	Gel	\$ 799.00
☐ Center Hole (CTHP)		1-Part	4" Opening x No Slot	☐ 801QC84	No Gel	\$ 642.00	☐ 801QC108	Gel	\$ 799.00
			5" Opening x No Slot	☐ 801QC85	No Gel	\$ 642.00	☐ 801QC109	Gel	\$ 799.00
			6" Opening x No Slot	☐ 801QC86	No Gel	\$ 642.00	☐ 801QC110	Gel	\$ 799.00

Back Cushion	\$
--------------	----

Laterals compatible with tall & narrow back only

☐ 801OD49	Padded Tall & Narrow Back Cushion - 9" Wide x 16" Tall	\$ 138.00
☐ 801OD50	Padded T-Back Cushion - 13" Wide x 11" Tall	\$ 138.00
☐ 801OD51	Padded Tall & Narrow Back Cushion - 11" Wide x 16" Tall	\$ 138.00
☐ 801OD52	Padded T-Back Cushion - 15" Wide x 11" Tall	\$ 138.00
☐ 801OD53	Padded Tall & Narrow Back Cushion - 13" Wide x 16" Tall	NC
☐ 801OD54	Padded T-Back Cushion - 17" Wide x 11" Tall	STD
☐ Custom Padded Back Cushion (select custom back cushion in next section)		BELOW

Custom Back Cushion	\$
---------------------	----

Laterals compatible with tall & narrow back only. 20" tall backs require backpost extenders.

☐ T-Back 	↳ 17" Width	☐ 801OD37	17" Wide x 14" Tall	☐ 801OD39	17" Wide x 18" Tall	\$ 642.00
		☐ 801OD38	17" Wide x 16" Tall	☐ 801OD40	17" Wide x 20" Tall	
☐ Tall & Narrow 	↳ 9" Width	☐ 801OD55	9" Wide x 18" Tall	☐ 801OD57	9" Wide x 20" Tall	\$ 642.00
	↳ 11" Width	☐ 801OD56	11" Wide x 18" Tall	☐ 801OD58	11" Wide x 20" Tall	
	↳ 13" Width	☐ 801OD45	13" Wide x 18" Tall	☐ 801OD46	13" Wide x 20" Tall	

☒ OPTION								RETAIL
Backposts								\$
☐ 801BP73	Single Bend Backposts (not compatible with cloth back) (decreases seat depth by ~2")						NC	
☐ 801BP74	Double Bend Backposts (not compatible with tall & narrow or cloth back)						NC	
☐ 801BP75	Triple Bend Backposts (compatible with cloth back only)						NC	
Leckey Chair Size								\$
Only select if ordering MC4000Leckey. Refer to Leckey size chart on last page. Leckey chair comes standard with Leckey chest and hip positioning belt.								
☐ 801SI30	Size 1 Leckey Chair	\$ 1,652.00	☐ 801SI32	Size 3 Leckey Chair	\$ 1,652.00			
☐ 801SI31	Size 2 Leckey Chair	\$ 1,652.00	☐ 801SI33	Size 4 Leckey Chair	\$ 1,652.00			
Footrest								
Footrests not compatible with MC4000Leckey. Seat mounted footrests compatible with MC4000TiltPed only.								
Footrest Style								\$
☐ 8013N02	Pediatric Adjustable Swing-Away Footrests (2) (not compatible with tilt models)						STD	
☐ Seat Mounted Footrests (2) (select sz)	☐ 801DI94	SM (12 - 14")	☐ 801DI95	MD (14 - 16")	☐ 801DI82	LG (17 - 20")	STD	
☐ 8014N07	Omit footrests						NC	
Footrest Accessories								\$
☐ 801DI74	Calf Strap	\$ 81.00	☐ 801DI78	Foot / Ankle Straps S (6.5 - 8")	\$ 545.00			
☐ 801DI75	BodyPoint Calf Panel M	\$ 255.00	☐ 801DI79	Foot / Ankle Straps M (7.5 - 9.25")	\$ 545.00			
☐ 801DI76	Heel Loops (2) (2" wide with mounting hardware)	\$ 131.00	☐ 801DI80	Foot / Ankle Straps L (9 - 11")	\$ 545.00			
☐ 801DI77	Foot / Ankle Straps XS (5.5 - 6.75")	\$ 545.00	☐ 801DI81	Foot / Ankle Straps XL (10.5 - 13")	\$ 545.00			
Armrest								
Armrests not compatible with MC4000Leckey. Swing-away locking side arms not compatible with MC4000TiltPed or side-mounted hip guides.								
Armrest Style								\$
☐ 801NP03	Cantilever Fold-Back Locking Arms						STD	
☐ 801NP05	Quick Release Cantilever Fold-Back Arms						\$ 109.00	
☐ 801NP07	Swing-Away Locking Side Arms						\$ 104.00	
Ampad Style								\$
Ampad selection for swing-away locking side arms not required. Arm troughs recommended with quick release arms due to width addition.								
☐ 801NP13	10" Arm Pads - Pair (compatible with swing-away locking side arms only)						NC	
☐ 801NP12	14" Arm Pads - Pair (not compatible with swing-away locking side arms)						STD	
☐ 801AM36	Arm Troughs - Pair (not compatible with swing-away locking side arms)						\$ 399.00	
☐ Arm Trough - Single (select side) (14" arm pad will be installed on side not selected)	☐ 801AM37	Left	☐ 801AM38	Right	\$ 200.00			
Headrest								\$
Headrest not compatible with MC4000Leckey								
☐ 801HE51	Adjustable Padded Headrest Assembly for Tilt Base with Padded Back Cushion						STD	
☐ 801HE50	Adjustable Padded Headrest Assembly for Non-Tilt Base with Padded Back Cushion						\$ 573.00	
Headrest Options								\$
☐ 801HE58	Quick Release Knobs for Headrest		☐ 801Y058	Large Headrest Pad (11.5 x 5.8")	\$ 145.00			
Upholstery Color								
Nuprodx Seat		\$			\$			
☐ 801UC40	Gray Color	STD	☐ 801UC59	Dark Blue Color	NC	☐ 801UC60	Camo Green	NC
Leckey Seat		\$			\$			
☐ 801UC61	Royal Blue	NC	☐ 801UC62	Pink				

✓ OPTION	DESCRIPTION					RETAIL
Accessories						
Non-Leckey			\$		\$	
Commode pan lowering brackets require commode pan. Quick release knobs require a selection of BodyPoint chest and/or hip positioning belt. Solid seat not compatible with splashguard.						
Chest harness compatible with tall & narrow back only. Laterals compatible with tall & narrow or cloth back only. Directional splashguard requires lowering brackets if ordering pan.						
☐ 801EJ12	Basic Chest Positioning Belt	\$	81.00	☐ 801GU8	Side-Mounted Hip Guide - Right (N/A with S/A arms)	\$ 273.00
☐ 801EJ09	BodyPoint Chest Positioning Belt M (24.5 - 34.5")	\$	275.00	☐ 801IG7	Under-Frame Splashguard (for 17" depth)	\$ 330.00
☐ 801EJ13	Basic Hip Positioning Seat Belt	\$	81.00	☐ 801IG8	Splashguard / Abductor (for front-facing slot only)	\$ 380.00
☐ 801EJ14	BodyPoint Hip Positioning Belt M (24.5 - 34.5")	\$	275.00	☐ 801KC18	Laterals - Pair	\$ 689.00
☐ 801BP72	Back Post Extenders (extends posts & back ht by 4")	\$	145.00	☐ 801KC31	Lateral - Left Side (1)	\$ 345.00
☐ 8013A12	Seat Depth Adapters (pair)	\$	159.00	☐ 801KC32	Lateral - Right Side (1)	\$ 345.00
☐ 8013A13	Quick Release Knobs for 1 BodyPoint Belt (1 pair)	\$	28.00	☐ 801IG9	Directional Under-Frame Splash Guard (for 17" width)	\$ 284.00
☐ 8013A22	Quick Release Knobs for 2 BodyPoint Belts (2 pair)	\$	55.00	☐ 801NH116	BodyPoint Chest Harness XS (9.5 - 11")	\$ 317.00
☐ 801PN16	Commode Pan with Rails (12" x 14" x 6")	\$	206.00	☐ 801NH117	BodyPoint Chest Harness S (11 - 13")	\$ 317.00
☐ 801PN17	Commode Pan Lowering Brackets (lowers pan by 4")	\$	131.00	☐ 801NH118	BodyPoint Chest Harness M (13 - 16")	\$ 317.00
☐ 801GU6	Side-Mounted Hip Guides - Pair (N/A with S/A arms)	\$	545.00	☐ 801NH119	BodyPoint Chest Harness L (16 - 19")	\$ 317.00
☐ 801GU7	Side-Mounted Hip Guide - Left (N/A with S/A arms)	\$	273.00	☐ 801NH120	BodyPoint Chest Harness XL (19 - 22")	\$ 317.00
MC4000Ped						\$
☐ Drop Rail (specify installation side) cushion may be rotated to match slot with drop side	☐ 8013A16	Left Side		☐ 8013A17	Right Side	
	☐ 8013A18	Front Side				
MC4020Ped						\$
☐ Drop Bar (specify installation side) cushion may be rotated to match slot with drop side	☐ 801SZ16	Left Side		☐ 801SZ17	Right Side	
	☐ 801SZ18	Front Side				
MC4000Leckey			\$		\$	

Color of accessories will match the chosen upholstery color

☐ 801A882	Leckey Abduction Belt (size 1 for size 1 or 2 Leckey chair, size 2 for size 3 or 4 Leckey chair)	\$	247.00
☐ 801HE57	Leckey Head Support	\$	261.00

Custom Requests

Please see the list below for common custom requests. Please note any of these in your order or contact NuprodX Custom Service.

Frame Size	Seat Cushions	Back Cushions	Seat-to-Floor Height	Footrest Angle / Length
Leg Height	Backpost Angle	Bridge Length	Armrest Pad Size	Hip Pad Size
Notes (optional)				

Leckey Size Chart

Seat width is 16" for all Leckey sizes

Product Sizing		Size 1		Size 2		Size 3		Size 4	
Age (Years)		1 - 5		4 - 9		8 - 14		12 - 18	
Advanced Bathchair		LAB/1		LAB/2		LAB/3		LAB/4	
		millimeters	inches	millimeters	inches	millimeters	inches	millimeters	inches
User Height	min	750	30	950	37	1150	45	1350	53
	max	1050	41	1250	49	1450	57	1650	65
Seat Depth		225	9	275	11	335	13	400	16
Back Support		350	14	400	16	460	18	550	22
Head Support		150	6	150	6	215	8.5	215	8.5
Calf Support		260	10	260	10	340	13	340	13

User Measurement Guide (optional)

If you would like us to have more information about your measurements, please fill out our user measurement guide below.

A	Chest / Lateral Width	
B	Hip Width	
C	Shoulder Height	
D	Seat Depth	
E	Lower Leg Length	

