



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		Province	
Buyer				PC		Phone	
Mark For				Email / Fax			

✓ OPTION	DESCRIPTION	TOTAL STANDARD PACKAGE	RETAIL
Base		\$	
☐ 801BS61	MC5000 (left-right sliding commode/shower/tub chair with crutch tips)	\$ 6,368.00	\$ 3,443.00
☐ 801BS62	MC5200 (heavy-duty left-right sliding commode/shower/tub chair with crutch tips)	\$ 6,644.00	\$ 3,443.00
☐ 801BS63	MC5100 (left-right sliding commode/shower/tub bench with crutch tips - 30" width - 2 piece)	\$ 3,181.00	\$ 3,181.00
☐ 801BS64	MC5300 (left-right sliding commode bench with crutch tips - custom width - 1 piece)	\$ 3,181.00	\$ 3,181.00

Base Options

MC5000 or MC5200 requires bridge, slide direction, and tub/shower base. MC5100 requires slide orientation. MC5300 not compatible with base options.

Bridge		\$	
☐ 801BV1	Bridge (connects the chair to the tub / shower base) (select slide orientation below)	\$ 1,435.00	

Slide Orientation

☐ 801EV03	Slide Orientation - Left (to user's left)	NC	☐ 801EV04	Slide Orientation - Right (to user's right)	NC
----------------------------------	---	----	----------------------------------	---	----

Tub / Shower Base

☐ 8014B1	Tub / Shower Base (connects to the chair via bridge)	\$ 1,490.00
☐ 8014B2	1.5" Offset Outer Legs	STD
☐ 8014B3	3" Offset Outer Legs (used when the difference between the frame width and tub base width is greater than 3")	\$ 83.00
☐ 8014B4	10.5" Inner Legs	STD
☐ 8014B5	7.5" Inner Legs (used when the difference between the inner and outer tub height is more than 4")	\$ 83.00

Base Type

Crutch Tips		\$	
☒ 8013A09	Crutch Tips (all four legs)	NC	

Frame Size

17" width compatible with non-heavy-duty base only. 20" width compatible with heavy-duty base only. MC5100 not compatible with 20x20 frame size. MC5300 compatible w/ all frame s

		\$	
☐ 801W70	17" W x 17" D Frame	STD	
☐ 801W71	17" W x 20" D Frame	\$ 414.00	
☐ 801W72	20" W x 17" D Frame (HD)	\$ 414.00	
☐ 801W74	20" W x 20" D Frame (HD)	\$ 276.00	

✓ OPTION	DESCRIPTION						RETAIL
Seat Cushion				\$		\$	
☐ 8" Wide Commode Opening with 4" Slot	☐ 801CU126	No Gel	STD	☐ 801CU125	Gel	\$	799.00
☐ 5.5" Wide Commode Opening with 4" Slot	☐ 801CU128	No Gel	\$ 117.00	☐ 801CU127	Gel	\$	799.00
☐ 5.5" Wide Commode Opening with No Slot	☐ 801CU130	No Gel	\$ 145.00	☐ 801CU129	Gel	\$	799.00
☐ No Commode Opening (solid seat)	☐ 801CU132	No Gel	\$ 164.00	☐ 801CU131	Gel	\$	799.00
☐ 801CU138	Custom Cushion (select custom seat cushion in next section)						BELOW

Custom Seat Cushion				\$		\$	
---------------------	--	--	--	----	--	----	--

Select style and size/gel option. Center hole with 8.0" opening compatible with 20x20 frame size only.

☐ Tear-Drop Open (TDO)		5" Opening x 3.0" Slot	☐ 801QC09	No Gel	\$ 642.00	☐ 801QC36	Gel	\$ 799.00
		6" Opening x 3.5" Slot	☐ 801QC10	No Gel	\$ 642.00	☐ 801QC37	Gel	\$ 799.00
		7" Opening x 4.0" Slot	☐ 801QC11	No Gel	\$ 642.00	☐ 801QC38	Gel	\$ 799.00
☐ Racetrack Open (RTO)		5" Opening x 3.0" Slot	☐ 801QC12	No Gel	\$ 642.00	☐ 801QC39	Gel	\$ 799.00
		6" Opening x 3.5" Slot	☐ 801QC13	No Gel	\$ 642.00	☐ 801QC40	Gel	\$ 799.00
		7" Opening x 4.0" Slot	☐ 801QC14	No Gel	\$ 642.00	☐ 801QC41	Gel	\$ 799.00
☐ Tear-Drop Closed (TDC)		5" Opening x No Slot	☐ 801QC15	No Gel	\$ 642.00	☐ 801QC42	Gel	\$ 799.00
		6" Opening x No Slot	☐ 801QC16	No Gel	\$ 642.00	☐ 801QC43	Gel	\$ 799.00
		7" Opening x No Slot	☐ 801QC17	No Gel	\$ 642.00	☐ 801QC44	Gel	\$ 799.00
☐ Racetrack Closed (RTC)		4" Opening x No Slot	☐ 801QC18	No Gel	\$ 642.00	☐ 801QC45	Gel	\$ 799.00
		5" Opening x No Slot	☐ 801QC19	No Gel	\$ 642.00	☐ 801QC46	Gel	\$ 799.00
		6" Opening x No Slot	☐ 801QC20	No Gel	\$ 642.00	☐ 801QC47	Gel	\$ 799.00
☐ Tear-Drop Side Cutout Left Side (TDSL)		7" Opening x No Slot	☐ 801QC21	No Gel	\$ 642.00	☐ 801QC48	Gel	\$ 799.00
		5" Opening x 3.0" Slot	☐ 801QC22	No Gel	\$ 642.00	☐ 801QC49	Gel	\$ 799.00
		6" Opening x 3.5" Slot	☐ 801QC23	No Gel	\$ 642.00	☐ 801QC50	Gel	\$ 799.00
☐ Tear-Drop Side Cutout Right Side (TDSR)		7" Opening x 4.0" Slot	☐ 801QC24	No Gel	\$ 642.00	☐ 801QC51	Gel	\$ 799.00
		5" Opening x 3.0" Slot	☐ 801QC25	No Gel	\$ 642.00	☐ 801QC52	Gel	\$ 799.00
		6" Opening x 3.5" Slot	☐ 801QC26	No Gel	\$ 642.00	☐ 801QC53	Gel	\$ 799.00
☐ Straight U-Slot (SUS)		7" Opening x 4.0" Slot	☐ 801QC27	No Gel	\$ 642.00	☐ 801QC54	Gel	\$ 799.00
		3" Opening x 3" Slot	☐ 801QC28	No Gel	\$ 642.00	☐ 801QC55	Gel	\$ 799.00
		4" Opening x 4" Slot	☐ 801QC29	No Gel	\$ 642.00	☐ 801QC56	Gel	\$ 799.00
☐ Center Hole (CTH)		5" Opening x 5" Slot	☐ 801QC30	No Gel	\$ 642.00	☐ 801QC57	Gel	\$ 799.00
		6" Opening x 6" Slot	☐ 801QC31	No Gel	\$ 642.00	☐ 801QC58	Gel	\$ 799.00
		5" Opening x No Slot	☐ 801QC32	No Gel	\$ 642.00	☐ 801QC59	Gel	\$ 799.00
		6" Opening x No Slot	☐ 801QC33	No Gel	\$ 642.00	☐ 801QC60	Gel	\$ 799.00
		7" Opening x No Slot	☐ 801QC34	No Gel	\$ 642.00	☐ 801QC61	Gel	\$ 799.00
		8" Opening x No Slot	☐ 801QC35	No Gel	\$ 642.00	☐ 801QC62	Gel	\$ 799.00

Seat Cushion Orientation							
--------------------------	--	--	--	--	--	--	--

Only available for cushions with a front-facing slot

☐ 801CU147	Slot Open Towards Front (if ordering a drop, slot side will match drop side)	STD
☐ 801CU148	Slot Open Towards Back (rotate cushion 180 degrees) (not compatible with drop)	NC

Back Cushion							
--------------	--	--	--	--	--	--	--

Laterals compatible with tall & narrow or cloth back only

17" Width Frame Size Only						\$	
☐ 801OD26	Padded Back Cushion - 17" Wide x 8" Tall (not compatible with headrest)						STD
☐ 801OD27	Padded T-Back Cushion - 17" Wide x 11" Tall						\$ 117.00
☐ 801OD29	Padded Tall & Narrow Back Cushion - 13" Wide x 16" Tall						NC

✓ OPTION	DESCRIPTION										RETAIL
20" Width Frame Size Only											\$
☐ 801OD30	Padded Back Cushion - 20" Wide x 8" Tall (not compatible with headrest)										STD
☐ 801OD31	Padded T-Back Cushion - 20" Wide x 11" Tall										\$ 117.00
☐ 801OD33	Padded Tall & Narrow Back Cushion - 16" Wide x 16" Tall										NC
Back Cushion (continued)											
Universal Frame Size											\$
☐ 801OD34	Cloth Back Assembly										NC
☐ Padded Custom Back Cushion (select custom back cushion in next section)											BELOW
Custom Back Cushion											\$
Select style and size based on frame size selected. T-back not compatible with laterals. 20" tall backs require backpost extenders.											
☐ T-Back		 17" Width	☐ 801OD37	14" Tall	☐ 801OD38	16" Tall	☐ 801OD39	18" Tall	☐ 801OD40	20" Tall	\$ 642.00
		 20" Width	☐ 801OD41	14" Tall	☐ 801OD42	16" Tall	☐ 801OD43	18" Tall	☐ 801OD44	20" Tall	
☐ Tall & Narrow		 17" Width	☐ 801OD45	13" Wide x 18" Tall			☐ 801OD46	13" Wide x 20" Tall			\$ 642.00
		 20" Width	☐ 801OD47	16" Wide x 18" Tall			☐ 801OD48	16" Wide x 20" Tall			
Backposts											\$
☐ 801BP73	Single Bend Backposts (not compatible with cloth back) (decreases seat depth by ~2")										NC
☐ 801BP74	Double Bend Backposts (not compatible with tall & narrow or cloth back)										NC
☐ 801BP75	Triple Bend Backposts (compatible with cloth back only)										NC
Bridge Length											\$
Only select bridge length if ordering the MC5000 or MC5200											
☐ 8012B1	13" Locking Bridge										STD
☐ 8012B5	9" Locking Bridge			\$ 192.00	☐ 8012B12	17" Locking Bridge			\$ 192.00		
☐ 8012B6	10" Locking Bridge			\$ 192.00	☐ 8012B13	18" Locking Bridge			\$ 192.00		
☐ 8012B7	11" Locking Bridge			\$ 192.00	☐ 8012B14	19" Locking Bridge			\$ 192.00		
☐ 8012B8	12" Locking Bridge			\$ 192.00	☐ 8012B2	20" Locking Bridge			\$ 192.00		
☐ 8012B9	14" Locking Bridge			\$ 192.00	☐ 8012B16	21" Locking Bridge			\$ 192.00		
☐ 8012B10	15" Locking Bridge			\$ 192.00	☐ 8012B17	22" Locking Bridge			\$ 192.00		
☐ 8012B11	16" Locking Bridge			\$ 192.00	☐ 8012B18	23" Locking Bridge			\$ 192.00		
☐ To Be Determined (provide "K" measurement in the bathroom fitting guide on last page)	Measure the distance from the center of the tub or shower drain / spout to center of where the customer desires roll-in section to be located (i.e. toilet).										TBD
Frame Width											
Only select frame width if ordering the MC5100 or MC5300											
MC5100 (2-Piece)											NC
☐ 801W93	30" Frame Width										
MC5300 (1-Piece)											NC
☐ 801W94	30" Frame Width	☐ 801W80	35" Frame Width	☐ 801W85	40" Frame Width	☐ 801W89	44" Frame Width				
☐ 801W76	31" Frame Width	☐ 801W81	36" Frame Width	☐ 801W86	41" Frame Width	☐ 801W90	45" Frame Width				
☐ 801W77	32" Frame Width	☐ 801W82	37" Frame Width	☐ 801W87	42" Frame Width	☐ 801W91	46" Frame Width				
☐ 801W78	33" Frame Width	☐ 801W83	38" Frame Width	☐ 801W88	43" Frame Width	☐ 801W92	47" Frame Width				
☐ 801W79	34" Frame Width	☐ 801W84	39" Frame Width								
Footrest											
Footrest Style											\$
Only one footrest style can be selected (adjustable swing-away or seat mounted).											
☐ 8014N04	Adjustable Swing-Away Footrests (2) (not compatible with MC5300)										\$ 366.00
☐ 801F139	Seat Mounted Footrests (2) - Small (12 - 14")										\$ 76.00
☐ 801F140	Seat Mounted Footrests (2) - Medium (14 - 16")										\$ 76.00
☐ 801F141	Seat Mounted Footrests (2) - Large (17 - 20")										\$ 76.00

✓ OPTION	DESCRIPTION					RETAIL
Footrest Accessories		\$				\$
☐ 801DI74	Calf Strap	\$	81.00	☐ 801DI78	Foot / Ankle Straps S (6.5 - 8")	\$ 545.00
☐ 801DI75	BodyPoint Calf Panel M (for 17" width only)	\$	255.00	☐ 801DI79	Foot / Ankle Straps M (7.5 - 9.25")	\$ 545.00
☐ 801DI83	BodyPoint Calf Panel L (for 20" width only)	\$	255.00	☐ 801DI80	Foot / Ankle Straps L (9 - 11")	\$ 545.00
☐ 801DI76	Heel Loops (2) (2" wide with mounting hardware)	\$	131.00	☐ 801DI81	Foot / Ankle Straps XL (10.5 - 13")	\$ 545.00
☐ 801DI77	Foot / Ankle Straps XS (5.5 - 6.75")	\$	545.00			

Armrest

Swing-away locking side arms not compatible with side-mounted hip guides

Armrest Style						\$
☐ 801NP03	Cantilever Fold-Back Locking Arms (17" depth frames only)					STD
☐ 801NP04	Long Cantilever Fold-Back Locking Arms (20" depth frames only)					STD
☐ 801NP05	Quick Release Cantilever Fold-Back Arms (17" depth frames only)					\$ 109.00
☐ 801NP09	Long Quick Release Cantilever Fold-Back Arms (20" depth frames only)					\$ 109.00
☐ 801NP07	Swing-Away Locking Side Arms (17" depth frames only)					\$ 104.00
☐ 801NP08	Long Swing-Away Locking Side Arms (20" depth frames only)					\$ 104.00

Armpad Style

Armpad selection for swing-away locking side arms not required. Arm troughs recommended with quick release arms due to width addition.

☐ 801NP13	10" Arm Pads - Pair (compatible with 17" swing-away locking side arms only)				NC
☐ 801NP12	14" Arm Pads - Pair (not compatible with 17" swing-away locking side arms)				STD
☐ 801AM36	Arm Troughs - Pair (not compatible with swing-away locking side arms)				\$ 399.00
☐ Arm Trough - Single (select side) (14" arm pad will be installed on side not selected)	☐ 801AM37	Left	☐ 801AM38	Right	\$ 200.00

Headrest

☐ 801HE50	Adjustable Padded Headrest Assembly for Padded Back Cushion					\$ 573.00
☐ 801HE61	Adjustable Padded Headrest Assembly for Cloth Back Cushion					\$ 731.00

Headrest Accessories

☐ 801HE58	Quick Release Knobs for Headrest	\$	40.00	☐ 801Y058	Large Headrest Pad (11.5 x 5.8")	\$ 145.00
----------------------------------	----------------------------------	----	-------	----------------------------------	----------------------------------	-----------

Upholstery Color

If ordering cloth back with camo green upholstery color, a gray color cloth back will be included

☐ 801UC40	Gray Color	STD	☐ 801UC59	Dark Blue Color	NC	☐ 801UC60	Camo Green	NC
----------------------------------	------------	-----	----------------------------------	-----------------	----	----------------------------------	------------	----

Accessories

Chest harness compatible with tall & narrow back only. Laterals compatible with tall & narrow or cloth back only. Directional splashguard requires lowering brackets if ordering pan.

Commode pan lowering brackets require commode pan. Quick release knobs require a selection of BodyPoint chest and/or hip positioning belt. Solid seat not compatible with splashgu:

☐ Drop Bar Rigidizer (specify installation side and frame size) cushion will be manufactured to match slot with drop side		☐ 8018A135	Left Side (for 17" depth)	☐ 801A137	Left Side (for 20" depth)	\$ 544.00
		☐ 8018A136	Right Side (for 17" depth)	☐ 801A138	Right Side (for 20" depth)	\$ 544.00
☐ Seat Cushion Overlay (n/a with tall & narrow back) (specify frame size)		☐ 801IU31	17x17 Frame	☐ 801IU33	20x20 Frame	\$ 421.00
		☐ 801IU32	17x20 Frame or 20x17 Frame			
☐ 8013A10	Side-Mounted Suction Cups (2) (for MC5100 & MC530)	\$ 192.00	☐ 801PN17	Commode Pan Lowering Brackets (lowers pan by 4")		\$ 131.00
☐ 8014N01	Side-Mounted Suction Cups - Additional Set (2) (for MC5000 & MC5200)	\$ 192.00	☐ 801GU6	Side-Mounted Hip Guides - Pair (N/A with S/A arms)		\$ 545.00
			☐ 801GU7	Side-Mounted Hip Guide - Left (N/A with S/A arms)		\$ 273.00
☐ 8013A11	Suction Cups / Crutch Tip Combo	\$ 255.00	☐ 801GU8	Side-Mounted Hip Guide - Right (N/A with S/A arms)		\$ 273.00
☐ 801EJ12	Basic Chest Positioning Belt	\$ 81.00	☐ 801IG7	Under-Frame Splashguard (for 17" depth)		\$ 330.00
☐ 801EJ09	BodyPoint Chest Positioning Belt M (24.5 - 34.5")	\$ 275.00	☐ 801IG10	Under-Frame Splashguard (for 20" depth)		\$ 330.00
☐ 801EJ10	BodyPoint Chest Positioning Belt L (31 - 41")	\$ 275.00	☐ 801IG8	Splashguard / Abductor (for front-facing slot only)		\$ 380.00
☐ 801EJ11	BodyPoint Chest Positioning Belt XL (32 - 52")	\$ 288.00	☐ 801KC18	Laterals - Pair		\$ 689.00
☐ 801EJ13	Basic Hip Positioning Seat Belt	\$ 81.00	☐ 801KC31	Lateral - Left Side (1)		\$ 345.00
☐ 801EJ14	BodyPoint Hip Positioning Belt M (24.5 - 34.5")	\$ 275.00	☐ 801KC32	Lateral - Right Side (1)		\$ 345.00
☐ 801EJ15	BodyPoint Hip Positioning Belt L (31 - 41")	\$ 275.00	☐ 801IG9	Directional Under-Frame Splashguard (for 17" width)		\$ 284.00
☐ 801EJ16	BodyPoint Hip Positioning Belt XL (32 - 52")	\$ 288.00	☐ 801IG11	Directional Under-Frame Splashguard (for 20" width)		\$ 284.00

✓ OPTION	DESCRIPTION		RETAIL
Accessories (Continued)		\$	\$
☐ 801BP72	Back Post Extenders (extends posts & back ht by 4")	\$ 145.00	☐ 801NH116 BodyPoint Chest Harness XS (9.5 - 11") \$ 317.00
☐ 8013A12	Seat Depth Adapters (pair)	\$ 159.00	☐ 801NH117 BodyPoint Chest Harness S (11 - 13") \$ 317.00
☐ 8013A13	Quick Release Knobs for 1 BodyPoint Belt (1 pair)	\$ 28.00	☐ 801NH118 BodyPoint Chest Harness M (13 - 16") \$ 317.00
☐ 8013A22	Quick Release Knobs for 2 BodyPoint Belts (2 pair)	\$ 55.00	☐ 801NH119 BodyPoint Chest Harness L (16 - 19") \$ 317.00
☐ 801PN16	Commode Pan with Rails (for 17" depth)	\$ 206.00	☐ 801NH120 BodyPoint Chest Harness XL (19 - 22") \$ 317.00
☐ 801PN18	Commode Pan with Rails (for 20" depth)	\$ 206.00	

Custom Requests

Please see the list below for common custom requests. Please note any of these in your order or contact NuprodX Custom Service.

Frame Size	Seat Cushions	Back Cushions	Seat-to-Floor Height	Footrest Angle / Length
Leg Height	Backpost Angle	Bridge Length	Armrest Pad Size	Hip Pad Size



Notes (optional)

User Measurement Guide & Bathroom Fitting Guide (optional)

If you would like us to have more information about your measurements, please fill out our user measurement guide below.

User Measurement Guide

A	Chest / Lateral Width	
B	Hip Width	
C	Shoulder Height	
D	Seat Depth	
E	Lower Leg Length	



Bathroom Fitting Guide

A	Overall length of bathtub		I	Distance from edge of bathtub ledge to closest obstacle (when toilet is NOT next to the tub)	
B	Length of bathtub floor		J	Distance from center of bathtub to edge of outside bathtub ledge (when toilet is NOT next to the tub)	
C	Width of bathtub floor (flat portion only)		K	Distance from center of toilet bowl to center of bathtub drain / spout	
D	Overall width of bathtub		L	Distance from top edge of bathtub to bathtub floor (inside tub height)	
E	Distance between edge of toilet bowl and bathtub		M	Distance from top edge of bathtub to bathroom floor (outside tub height)	
F	Widest portion of toilet bowl (including bidet with if one is present)		N	Distance from top of toilet bowl to floor (with seat / lid in raised position)	
G	Distance between side edge of toilet bowl and closest object / wall		O	Doorway Width	
H	Distance from front edge of toilet bowl to closest object / wall				

If ordering for a shower, please include the following measurements

1	Height of the shower threshold on outside of shower		4	Overall width of the shower (left to right)	
2	Width of the shower entry		5	Overall depth of the shower (front to back)	
3	Height of the shower threshold on the inside of shower				



Customer				Ship To			
Account			Name				
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer			ZIP		Phone		
Mark For			Email / Fax				

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

✓ OPTION	HCPCS	DESCRIPTION	STANDARD CONFIGURATION	RETAIL		
Base				\$		
☐ 801BS61	E1399/E0240	MC5000 (left-right sliding commode/shower/tub chair with crutch tips)	\$ 4,615.00	\$ 2,495.00		
☐ 801BS62	E1399/E0240	MC5200 (heavy-duty left-right sliding commode/shower/tub chair with crutch tips)	\$ 4,815.00	\$ 2,495.00		
☐ 801BS63	E1399/E0240	MC5100 (left-right sliding commode/shower/tub bench with crutch tips - 30" width - 2 piece)	\$ 2,305.00	\$ 2,305.00		
☐ 801BS64	E1399/E0240	MC5300 (left-right sliding commode bench with crutch tips - custom width - 1 piece)	\$ 2,305.00	\$ 2,305.00		
Base Options						
MC5000 or MC5200 requires bridge, slide direction, and tub/shower base. MC5100 requires slide orientation. MC5300 not compatible with base options.						
Tilt				\$		
☒ 801TT14	Omit Tilt on Base			NC		
Bridge				\$		
☐ 801BV1	E1399	Bridge (connects the chair to the tub / shower base) (select slide orientation below)		\$ 1,040.00		
Slide Orientation				\$		
☐ 801EV03	Slide Orientation - Left (to user's left)		NC	☐ 801EV04	Slide Orientation - Right (to user's right)	NC
Tub / Shower Base				\$		
☐ 8014B1	E1399	Tub / Shower Base (connects to the chair via bridge)		\$ 1,080.00		
☐ NEW	1.5" Offset Outer Legs (for standard tubs)			STD		
☐ NEW	3" Offset Outer Legs (for narrow/wide tubs)			\$ 60.00		
☐ NEW	10.5" Inner Legs			STD		
☐ NEW	7.5" Inner Legs (inner tub ht 4" or more < outer)			\$ 60.00		
Base Type						
Crutch Tips				\$		
☒ 8013A09	Crutch Tips (all four legs)			NC		
Frame Size					\$	
17" width compatible with non-heavy-duty base only. 20" width compatible with heavy-duty base only. MC5100 not compatible with 20x20 frame size. MC5300 compatible w/ all frame s						
☐ 801W70	17" W x 17" D Frame and Cushion			STD		
☐ 801W71	17" W x 20" D Frame and Cushion			\$ 300.00		
☐ 801W72	20" W x 17" D Frame and Cushion (HD)			\$ 300.00		
☐ 801W74	20" W x 20" D Frame and Cushion (HD)			\$ 200.00		

Select style and size/gel option. Center hole with 8.0" opening compatible with 20x20 frame size only.

Only available for cushions with a front-facing slot

Laterals compatible with tall & narrow or cloth back only

17" Width Frame Size Only		\$
☐ 8010D26	Padded Back Cushion - 17" Wide x 8" Tall (not compatible with headrest)	STD
☐ 8010D27	Padded T-Back Cushion - 17" Wide x 11" Tall	\$ 85.00
☐ 8010D29	Padded Tall & Narrow Back Cushion - 13" Wide x 16" Tall	NC
20" Width Frame Size Only		\$
☐ 8010D30	Padded Back Cushion - 20" Wide x 8" Tall (not compatible with headrest)	STD

☐ 801OD31	Padded T-Back Cushion - 20" Wide x 11" Tall	\$ 85.00
☐ 801OD33	Padded Tall & Narrow Back Cushion - 16" Wide x 16" Tall	NC

MK-130434 Rev B © 12.2025

Page 2 of 5

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Back Cushion (continued)

Universal Frame Size	\$
----------------------	----

☐ 801OD34	Cloth Back Assembly	NC
☐ Padded Custom Back Cushion (select custom back cushion in next section)		BELOW

Custom Back Cushion	\$
---------------------	----

Select style and size based on frame size selected. T-back not compatible with laterals. 20" tall backs require backpost extenders.

☐ T-Back 	☐ 17" Width ☐ 801OD37 14" Tall ☐ 801OD38 16" Tall ☐ 801OD39 18" Tall ☐ 801OD40 20" Tall	\$ 465.00
	☐ 20" Width ☐ 801OD41 14" Tall ☐ 801OD42 16" Tall ☐ 801OD43 18" Tall ☐ 801OD44 20" Tall	
☐ Tall & Narrow 	☐ 17" Width ☐ 801OD45 13" Wide x 18" Tall ☐ 801OD46 13" Wide x 20" Tall	\$ 465.00
	☐ 20" Width ☐ 801OD47 16" Wide x 18" Tall ☐ 801OD48 16" Wide x 20" Tall	

Backposts	\$
-----------	----

☐ NEW	Single Bend Backposts (not compatible with cloth back) (decreases seat depth by ~2")	NC
☐ NEW	Double Bend Backposts (not compatible with tall & narrow or cloth back)	NC
☐ NEW	Triple Bend Backposts (compatible with cloth back only)	NC

Bridge Length	\$	\$
---------------	----	----

Only select bridge length if ordering the MC5000 or MC5200

☐ 8012B1	13" Locking Bridge		STD
☐ 8012B5	9" Locking Bridge	\$ 139.00	☐ 8012B12 17" Locking Bridge \$ 139.00
☐ 8012B6	10" Locking Bridge	\$ 139.00	☐ 8012B13 18" Locking Bridge \$ 139.00
☐ 8012B7	11" Locking Bridge	\$ 139.00	☐ 8012B14 19" Locking Bridge \$ 139.00
☐ 8012B8	12" Locking Bridge	\$ 139.00	☐ 8012B2 20" Locking Bridge \$ 139.00
☐ 8012B9	14" Locking Bridge	\$ 139.00	☐ 8012B16 21" Locking Bridge \$ 139.00
☐ 8012B10	15" Locking Bridge	\$ 139.00	☐ 8012B17 22" Locking Bridge \$ 139.00
☐ 8012B11	16" Locking Bridge	\$ 139.00	☐ 8012B18 23" Locking Bridge \$ 139.00

☐ To Be Determined (provide "K" measurement in the bathroom fitting guide on last page)		TBD
Measure the distance from the center of the tub or shower drain / spout to center of where the customer desires roll-in section to be located (i.e. toilet).		

Frame Width

Only select frame width if ordering the MC5100 or MC5300

MC5100 (2-Piece)	NC
------------------	----

☐ 801W93	30" Frame Width
---------------------------------	-----------------

MC5300 (1-Piece)	NC
------------------	----

☐ 801W94	30" Frame Width	☐ 801W80	35" Frame Width	☐ 801W85	40" Frame Width	☐ 801W89	44" Frame Width
☐ 801W76	31" Frame Width	☐ 801W81	36" Frame Width	☐ 801W86	41" Frame Width	☐ 801W90	45" Frame Width
☐ 801W77	32" Frame Width	☐ 801W82	37" Frame Width	☐ 801W87	42" Frame Width	☐ 801W91	46" Frame Width
☐ 801W78	33" Frame Width	☐ 801W83	38" Frame Width	☐ 801W88	43" Frame Width	☐ 801W92	47" Frame Width
☐ 801W79	34" Frame Width	☐ 801W84	39" Frame Width				

Footrest

Footrest Style	\$
----------------	----

Only one footrest style can be selected (adjustable swing-away or seat mounted).

☐ 8014N04	Adjustable Swing-Away Footrests (2) (not compatible with MC5300)	\$ 265.00
☐ 801F139	Seat Mounted Footrests (2) - Small (12 - 14")	\$ 55.00
☐ 801F140	Seat Mounted Footrests (2) - Medium (14 - 16")	\$ 55.00
☐ 801F141	Seat Mounted Footrests (2) - Large (17 - 20")	\$ 55.00

Footrest Accessories	\$	\$
----------------------	----	----

☐ 801DI74	Calf Strap	\$ 59.00	☐ 801DI78	Foot / Ankle Straps S (6.5 - 8")	\$ 395.00
----------------------------------	------------	----------	----------------------------------	----------------------------------	-----------

☐ 801DI75	BodyPoint Calf Panel M (for 17" width only)	\$ 185.00	☐ 801DI79	Foot / Ankle Straps M (7.5 - 9.25")	\$ 395.00
☐ 801DI83	BodyPoint Calf Panel L (for 20" width only)	\$ 185.00	☐ 801DI80	Foot / Ankle Straps L (9 - 11")	\$ 395.00
☐ 801DI76	Heel Loops (2) (2" wide with mounting hardware)	\$ 95.00	☐ 801DI81	Foot / Ankle Straps XL (10.5 - 13")	\$ 395.00
☐ 801DI77	Foot / Ankle Straps XS (5.5 - 6.75")	\$ 395.00			

MK-130434 Rev B © 12.2025

Page 3 of 5

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Armrest

Swing-away locking side arms not compatible with side-mounted hip guides

Armrest Style		\$
☐ 801NP03	Cantilever Fold-Back Locking Arms (17" depth frames only)	STD
☐ 801NP04	Long Cantilever Fold-Back Locking Arms (20" depth frames only)	STD
☐ 801NP05	Quick Release Cantilever Fold-Back Arms (17" depth frames only)	\$ 79.00
☐ 801NP09	Long Quick Release Cantilever Fold-Back Arms (20" depth frames only)	\$ 79.00
☐ 801NP07	Swing-Away Locking Side Arms (17" depth frames only)	\$ 75.00
☐ 801NP08	Long Swing-Away Locking Side Arms (20" depth frames only)	\$ 75.00

Armpad Style	\$
--------------	----

Armpad selection for swing-away locking side arms not required. Arm troughs recommended with quick release arms due to width addition.

☐ 801NP13	10" Arm Pads - Pair (compatible with 17" swing-away locking side arms only)				NC
☐ 801NP12	14" Arm Pads - Pair (not compatible with 17" swing-away locking side arms)				STD
☐ 801AM36	Arm Troughs - Pair (not compatible with swing-away locking side arms)				\$ 290.00
☐ Arm Trough - Single (select side) (14" arm pad will be installed on side not selected)	☐ 801AM37	Left	☐ 801AM38	Right	\$ 145.00

Headrest	\$
----------	----

☐ 801HE50	Adjustable Padded Headrest Assembly for Padded Back Cushion	\$ 415.00
☐ 801HE61	Adjustable Padded Headrest Assembly for Cloth Back Cushion	\$ 530.00

Headrest Accessories	\$
----------------------	----

☐ 801HE58	Quick Release Knobs for Headrest	\$ 29.00	☐ 801Y058	Large Headrest Pad (11.5 x 5.8")	\$ 105.00
----------------------------------	----------------------------------	----------	----------------------------------	----------------------------------	-----------

Upholstery Color	\$	\$	\$
------------------	----	----	----

If ordering cloth back with camo green upholstery color, a gray color cloth back will be included

☐ 801UC40	Gray Color	STD	☐ 801UC59	Dark Blue Color	NC	☐ 801UC60	Camo Green	NC
----------------------------------	------------	-----	----------------------------------	-----------------	----	----------------------------------	------------	----

Accessories	\$	\$
-------------	----	----

Chest harness compatible with tall & narrow back only. Laterals compatible with tall & narrow or cloth back only. Directional splashguard requires lowering brackets if ordering pan.

Commode pan lowering brackets require commode pan. Quick release knobs require a selection of BodyPoint chest and/or hip positioning belt. Solid seat not compatible with splashguard

☐ Drop Bar Rigidizer (specify installation side and frame size) cushion will be manufactured to match slot with drop side	☐ 8018A135	Left Side (for 17" depth)	☐ 801A137	Left Side (for 20" depth)	\$ 394.00
	☐ 8018A136	Right Side (for 17" depth)	☐ 801A138	Right Side (for 20" depth)	\$ 394.00
☐ Seat Cushion Overlay (n/a with tall & narrow back) (specify frame size)	☐ 801IU31	17x17 Frame	☐ 801IU33	20x20 Frame	\$ 305.00
	☐ 801IU32	17x20 Frame or 20x17 Frame			
☐ 8013A10	Side-Mounted Suction Cups (2) (for MC5100 & MC530)	\$ 139.00	☐ 801PN17	Commode Pan Lowering Brackets (lowers pan by 4")	\$ 95.00
☐ 8014N01	Side-Mounted Suction Cups - Additional Set (2) (for MC5000 & MC5200)	\$ 139.00	☐ 801GU6	Side-Mounted Hip Guides - Pair (N/A with S/A arms)	\$ 396.00
			☐ 801GU7	Side-Mounted Hip Guide - Left (N/A with S/A arms)	\$ 198.00
☐ 8013A11	Suction Cups / Crutch Tip Combo	\$ 185.00	☐ 801GU8	Side-Mounted Hip Guide - Right (N/A with S/A arms)	\$ 198.00
☐ 801EJ12	Basic Chest Positioning Belt	\$ 59.00	☐ 801IG7	Under-Frame Splashguard (for 17" depth)	\$ 239.00
☐ 801EJ09	BodyPoint Chest Positioning Belt M (24.5 - 34.5")	\$ 199.00	☐ 801IG10	Under-Frame Splashguard (for 20" depth)	\$ 239.00
☐ 801EJ10	BodyPoint Chest Positioning Belt L (31 - 41")	\$ 199.00	☐ 801IG8	Splashguard / Abductor (for front-facing slot only)	\$ 275.00
☐ 801EJ11	BodyPoint Chest Positioning Belt XL (32 - 52")	\$ 209.00	☐ 801KC18	Laterals - Pair	\$ 500.00
☐ 801EJ13	Basic Hip Positioning Seat Belt	\$ 59.00	☐ 801KC31	Lateral - Left Side (1)	\$ 250.00
☐ 801EJ14	BodyPoint Hip Positioning Belt M (24.5 - 34.5")	\$ 199.00	☐ 801KC32	Lateral - Right Side (1)	\$ 250.00
☐ 801EJ15	BodyPoint Hip Positioning Belt L (31 - 41")	\$ 199.00	☐ 801IG9	Directional Under-Frame Splashguard (for 17" width)	\$ 206.00
☐ 801EJ16	BodyPoint Hip Positioning Belt XL (32 - 52")	\$ 209.00	☐ 801IG11	Directional Under-Frame Splashguard (for 20" width)	\$ 206.00
☐ 801BP72	Back Post Extenders (extends posts & back ht by 4")	\$ 105.00	☐ 801NH116	BodyPoint Chest Harness XS (9.5 - 11")	\$ 230.00

☐ 8013A12	Seat Depth Adapters (pair)	\$ 115.00	☐ 801NH117	BodyPoint Chest Harness S (11 - 13")	\$ 230.00
☐ 8013A13	Quick Release Knobs for 1 BodyPoint Belt (1 pair)	\$ 20.00	☐ 801NH118	BodyPoint Chest Harness M (13 - 16")	\$ 230.00
☐ 8013A22	Quick Release Knobs for 2 BodyPoint Belts (2 pair)	\$ 40.00	☐ 801NH119	BodyPoint Chest Harness L (16 - 19")	\$ 230.00
☐ 801PN16	Commode Pan with Rails (for 17" depth)	\$ 149.00	☐ 801NH120	BodyPoint Chest Harness XL (19 - 22")	\$ 230.00
☐ 801PN18	Commode Pan with Rails (for 20" depth)	\$ 149.00			

MK-130434 Rev B © 12.2025

Page 4 of 5

Custom Requests

Please see the list below for common custom requests. Please note any of these in your order or contact NuprodX Custom Service.

Frame Size	Seat Cushions	Back Cushions	Seat-to-Floor Height	Footrest Angle / Length
Leg Height	Backpost Angle	Bridge Length	Armrest Pad Size	Hip Pad Size



Notes (optional)

User Measurement Guide & Bathroom Fitting Guide (optional)

If you would like us to have more information about your measurements, please fill out our user measurement guide below.

User Measurement Guide

A	Chest / Lateral Width	
B	Hip Width	
C	Shoulder Height	
D	Seat Depth	
E	Lower Leg Length	
F	Seat to Floor Height	



Bathroom Fitting Guide

A	Overall length of bathtub		I	Distance from edge of bathtub ledge to closest obstacle (when toilet is NOT next to the tub)	
B	Length of bathtub floor		J	Distance from center of bathtub to edge of outside bathtub ledge (when toilet is NOT next to the tub)	
C	Width of bathtub floor (flat portion only)		K	Distance from center of toilet bowl to center of bathtub drain / spout	
D	Overall width of bathtub		L	Distance from top edge of bathtub to bathtub floor (inside tub height)	
E	Distance between edge of toilet bowl and bathtub				
F	Widest portion of toilet bowl (including bidet with if one is present)				

G	Distance between side edge of toilet bowl and closest object / wall		M	Distance from top edge of bathtub to bathroom floor (outside tub height)	
H	Distance from front edge of toilet bowl to closest object / wall		N	Distance from top of toilet bowl to floor (with seat / lid in raised position)	
			O	Doorway Width	

If ordering for a shower, please include the following measurements

1	Height of the shower threshold on outside of shower		4	Overall width of the shower (left to right)	
2	Width of the shower entry		5	Overall depth of the shower (front to back)	
3	Height of the shower threshold on the inside of shower				

Customer Service Tel: 707.934.8266 | Email: NuprodxCustomerService@sunmed.com | Fax: 707.934.8269 | Web: www.SunriseMedical.com

MK-130434 Rev B © 12.2025

