



Serial No:



ORDER FORM

Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg User height
Build4me			

Product Style

ELAXZ010	Elite	\$5,321.00
ELAXA062	Standard Frame	Std
ELAXA063	Standard Frame with 25mm seat rail	\$76.00
ELAXA280	Heavy Duty Frame (Spinergy R10 Wheels also required)	\$607.00
ELAXA145	Rugby League Upgrade	\$481.00

Max User Weight 125kg (275lbs)

Max User Weight ELAXA280 Heavy duty frame 158kg (348lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD Approval

XW009	Quote only (No CAD For Approval)	
XW010	Order CAD Drawing For Approval (Includes 2 CAD Amendments)	\$300.00 FLAT FEE
XW011	Order (No CAD Required)	
XW012	Quote & CAD For Approval (Includes 2 CAD Amendments)	\$550.00 FLAT FEE

Key Measurements

XY010	Seat Width (Inside to inside of sideguards)	
	mm (Max - 500mm)	
XY020	Seat Length	
	mm (Max - 500mm)	
XY030	Extra Seat Length (Provides additional leg support & hand grip)	
	mm	
XY050	Front frame height	(Min - 360mm)
	mm	(Max - 630mm)
XY060	Rear frame height	(Min - 360mm)
	mm	(Max - 630mm)
XY070	Frame Length	
	mm	(Min = XY020 + XY030 + 140mm)
XA020	Ergonomic seat	\$92.00
	Specify Length	
	mm	(Standard = 120mm)

Front Frame and Seat

Front Frame Options

XY080	Seat taper design	XA130	V-Seat design	\$245.00
Specify taper width		Specify taper width		
mm		mm		
		NOTE: Inside to inside wheel (XY470) required		

Advanced Front Frame Options

XY082	Seat taper position			\$91.00
Specify measurement				
mm				
XA131	V-seat and taper design			\$331.00
Specify inside width at front of sideguard (A)	Specify taper width (B)			
mm	mm			
	NOTE: Inside to inside wheel (XY470) required			
XA136	Parallel taper design (Not with seat taper XY080)			\$91.00
Specify taper width				
mm				
XA135	Asymmetrical taper design			\$387.00
Specify taper width (A)	Specify offset (B) (inside of sideguard to inside of taper width)			
mm	mm			
Specify offset taper side				
Left Right				

Seat Options

XR040	Aluminium seat plate (riveted to frame)			\$153.00
Specify length of plate	(Std XY020 - 20mm)			
mm				
XR050	Carbon fibre seat plate (riveted to frame)			\$391.00
Specify length of plate	(Std XY020 - 20mm)			
mm				
XA260	Drop seat design (Only with either XR040 or 050)			\$213.00
Specify drop				
mm	(Std = 20mm)			

Footrest

ELAXB080	Fixed position footrest	Std
ELAXB130	Clamp on adjustable position footrest (100mm range)	\$233.00
ELAXB085	Adjustable position footrest (50mm range)	\$142.00
	Front footrest height	
	mm (min = 35mm fixed / 40mm adjustable)	
	Rear footrest height	
	mm (min = 35mm fixed / 40mm adjustable)	
	Footrest position	
	(front of footplate to castor pot)	
	Footrest width (Not ABS)	
	mm (Std= 200mm)	
	Footrest length (Not ABS)	
	mm (Std= 250mm)	

Fixed Footrest Plate Options

XB180	Aluminium footplate (2.5mm)	Std
XB190	Carbon fibre footplate	\$254.00
XB100	ABS Moulded foot tray	\$135.00
	Specify length	
	mm Max = 380mm (Full) Std = 330mm	
ELAXB140	Spacer kit	\$21.00
	(Allows for 25mm increase in height)	
ELAXB145	Multi position footrest tubes	\$149.00

Adjustable Footrest Plate Options

XB180	Aluminium footplate (2.5mm)	No cost
XB100	ABS moulded foot tray	\$149.00
	Specify length	
	mm Max = 380mm (Full) Std = 330mm	

Additional Footrest Plate Options

XB240	Aluminium raised edges	\$87.00
	(Cost per edge)	
	Front Specify height	Specify Depth
	Back (Max = 60mm)	
	Left	
	Right	
XB250	Carbon Fibre raised edge	\$336.00
	(1 side only)	
	Front Specify height	
	Back (Max = 45mm)	

Castor pin width

XY330	Castor pin width
	mm MAX 320mm
	Optimise measurements

Camber Bar, Wing and Anti-Tip

ELAXY510	<p>Fixed centre of gravity</p> <p>Specify measurement</p> <p style="text-align: center;">mm</p>	<p>(Note: Dimension is taken from front of backrest to centre of axle)</p>	No cost
----------	--	--	----------------

ELAXY515	<p>Backrest COG adjustment</p> <p>(not with welded sides) (+/- 20mm)</p> <p>Specify measurement</p> <p style="text-align: center;">mm</p>	\$182.00
----------	--	-----------------

Specify Camber Angle

XY390	<p>Camber angle</p> <p>16 18 20 22</p>
-------	---

Wing Options

ELAXA110	<p>Round design basketball wing</p>	No cost
ELAXA120	<p>Straight sided design basketball wing</p>	\$193.00
ELAXA125	<p>Arrow design basketball wing</p>	\$193.00

Anti-Tip options

ELAXJ040	<p>Single fixed anti-tip</p>	No cost
ELAXJ060	<p>Double fixed anti-tip</p> <p>(With high mount joining bar)</p> <p>Specify castor pin width</p> <p style="text-align: center;">mm</p> <p>Optimise measurements (Wide as possible - std if not indicated)</p>	\$274.00

Sideguards

ELAXE230	Welded aluminium sport (Not with backrest COG Adj - ELAXY510)	\$413.00
ELAXE250	Welded aluminium tubes only (Not with backrest COG Adj - ELAXY510)	\$193.00
XE310	High design upgrade (Only with option ELAXE230)	\$194.00
XE030	Screw on aluminium	\$204.00
XE050	Screw on carbon fibre (Not recommended for basketball use)	\$456.00

Sideguard Design (Choose one)

Straight design

Round design

Follow tyre

Sideguard Measurement

Sideguard length

mm

Follow tyre design

25mm above tyre

Sideguard height (back)

mm

Individual above tyre

mm

Sideguard height (front)

(Optional for max height players)

mm

Level with tyre

Backrest Height

XY230

Fixed backrest height

Specify height

mm

XY245

Adjustable backrest height

(+/-25mm)

Specify height

mm

\$153.00

Ergonomic Options

XY360 Backrest angle

-5 -3 0 +3 +5

Individual

Deg

(Max -10/+10 Deg)

XG100

Ergonomic backrest

Vertical height

mm

Rearward position

mm

\$92.00

XG110 Inward taper

(not with ELAXY245)

Specify width or per side

mm

\$92.00

XG120 Outward flare

(not with ELAXY245)

Specify width or per side

mm

\$92.00

Castor Options

XC160	Sport 61mm Castor (Black)	XC181	Bones Reds bearings (All castor bearings)	
XC200	Sport micro fork	XC182	Bones Swiss Ceramic bearings (All castor bearings)	\$322.00

Rear Wheels

Rear Wheel Size

XY420	24" Wheels	XY430	25" Wheels	XY440	26" Wheels	XY450	700C Wheels
-------	------------	-------	------------	-------	------------	-------	-------------

Rear Wheel Gap

XY460	Rear Wheel Gap	XY470	Inside to inside of wheel
15	20	25	30
Individual	mm	mm	

NOTE: Required for V seat & taper

Rear Wheel Options

XK030	Equalizer 36/48 Spoke - No cost	XK330	Spinergy SLX R10 heavy duty 24 spoke (18 degree camber only)	\$1,527.00
XK080	Spinergy SLX 24 Spoke Spinergy			\$867.00
XK101	XSLX 24 Spoke + Caps			\$914.00
XK040	Supply without wheels and axle pins	Wheel type:	Tyre type:	

Spinergy Colour Options

Black Spokes	White Spokes	Yellow Spokes	Red Spokes	Green Spokes
Blue Spokes	Pink Spokes	Purple Spokes	Orange Spokes	
XK180	Individual spoke design (Specify colour)	Inside: Outside:		\$203.00

Spinergy Anodise Options

XK190	Black hub (Std option)	XK200	Silver hub	\$106.00
XK201	Red hub/rims			\$316.00
XK202	Purple hub/rims			\$316.00
XK203	Blue hub/rims			\$316.00
XK204	Gold hub/rims			\$316.00

Spare Wheel

ELAXK500	Spare wheel spinergy (Built up + axle pin)	\$978.00
XK520	Spare wheel upcharge - Titanium rims	\$227.00
XK522	Spare wheel upcharge - Remove excess tabs	\$60.00
XK526	Spare wheel upcharge - R10 Hub	\$350.00

Handrim Position/Options

XY480	Narrow set handrim	
XY485	Wide set handrim	
XL090	Remove excess mounting tabs (recommended for small wheel gaps)	\$117.00
XL080	Super narrow set handrim (only with XL010 & XL070 and includes XL090)	\$235.00

Handrim Style

XL010	Aluminium hard anodised 6 tab	Std
XL040	Titanium 6 tab (Not with 700C)	\$278.00
XL070	Aluminium hard anodised 12 tab (24" not available)	\$117.00
XL050	Stainless steel 6 tab	\$146.00

Tyre Options

XM070	Kenda Kaliente	Red	Blue	Std
XM150	IRC Exeracer Pro	(Red only)		\$96.00
XM080	Tufo MS3 tubeless	(Blue only)		\$96.00
XM095	Kenda Wide	(Grey only)		\$41.00
XM096	Schwable Dunk King			\$243.00
	(Not available in 700c)			

Inner Tube Valve Options

XY500	Schraeder inner tube valve	No Cost
XY490	Presta inner tube valve (Standard if no option selected)	Std

Finish Options

XN010 Powder coated frame

No Cost

Frame

Gloss

- Black gloss 103
- White gloss 100
- Grey gloss 132
- Cream gloss 135
- Red gloss 108
- Yellow gloss 120
- Forest green gloss 143
- Green gloss 107
- Dark blue gloss 151
- Royal blue gloss 154
- Blue gloss 174
- Duckegg Blue gloss 141
- Pink gloss 188
- Orange gloss 180
- Purple gloss 157

Frame

Matt

- Black matt 105
- White matt 101
- Grey matt 131
- Light grey matt 133
- Navy blue matt 152

Frame

Special

- Chrome mirror 300
- Sparkle silver 170
- Metallic black matt 104
- Metallic grey matt 130
- Black sandpaper 201
- Grey sandpaper 202
- White sandpaper 200

XN150 Powder coated components

(Frame colour if not specified)

No Cost

Specify component colour code from above list

ELAXN030

Brushed Frame

(All parts black unless XN160 chosen)

\$383.00

XN160

Brushed Components

\$145.00

ELAXN050

Polished Frame

(All parts black unless XN170 chosen)

\$1,013.00

XN170

Polished Components

\$213.00

Accent Colour

Note: Accent colour includes forks only

Black accent

Silver accent

Red accent

Blue accent

Purple accent

Upholstery

Seat upholstery

XP040 Seat sling strap adjustable (Std) No seat sling (Seat Plate) No Cost

Backrest upholstery

XP055 Airtech sport backrest strap adjustable (Black) XP070 Backrest upgrade Leather Corners \$132.00

Backrest Stitching options

Silver	Black	Blue	Individual:
Red	Orange	Purple	

RGK Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

Model Logo colour

White (std)	Green	Blue	Individual:
Red	Orange	Purple	

XP170 Individual backrest embroidery \$150.00
(PLEASE PROVIDE LOGO AT TIME OF ORDER)

Cushion

XP080 Cushion Foam \$92.00
(Select density and size below)

Soft (CMS040)		Medium (CMS050)		
50mm	75mm	100mm	Other	mm

XP340 Hard Plastazote cushion upgrade \$86.00
(Indicate max height)

580mm 630mm

Cushion Cover Options

XP085	Nylon cushion cover (Black)	Std
XP380	Neoprene cushion cover (Black)	\$98.00
XP090	Air Lite cushion cover (Black)	\$71.00

Straps

XS050	100mm Pro adjustable velcro strap	\$131.00
XS055	50mm Pro adjustable velcro strap	\$131.00
XS057	50mm Pro knee strap (measurement from tube to tube over knees)	\$131.00

Strap 1	Type	Measurement		
Strap 2	Type	Measurement	Specify strap length	Distance from back post
Strap 3 (Additional)	Type	Measurement	mm	mm

XS020	Sport ratchet strap (Click strap) (11" x 2" pad)	Quantity:	1	2	\$89.00
-------	---	-----------	---	---	---------

XS021	Sport ratchet strap (Single)	Quantity:	1	2	\$90.00
-------	------------------------------	-----------	---	---	---------

XS025	Pelvi Loc 2DA ratchet strap (For thighs & knees)	Quantity:	1	2	\$233.00
	Select size:				
	Medium (9")	Large (11")			

XS030	Pelvi Loc 3DA ratchet strap (For hips)	Quantity:	1	2	\$250.00
	Select size:				
	XS (8")	S (12")	M (14")	L (17")	

XS035	M2 Ratchet strap (For thighs and knees)	Quantity:	1	2	\$149.00
	6" 8" 10" 12"				

XS036	M2 Ratchet strap (Single) (For Feet)	Quantity:	1	2	\$149.00
	6" 8" 10" 12"				

Strap 1	Type	Strap 2	Type	Strap 3	Type
				Specify distance from back post	
				mm	

	Specify distance from back post	Specify distance from back post		
	mm	mm		

XS040	Fixed mounting points (One set per click strap)	x1	x2	x3	\$93.00
-------	--	----	----	----	---------

XS045	Clamped mounting points (One set per click strap)	x1	x2	x3	\$96.00
-------	--	----	----	----	---------

Frame options

XR060	Calf strap mounting loops	\$154.00
-------	---------------------------	----------

XR065	Strap location loops	\$111.00
-------	----------------------	----------

Accessories

XR030	Triple wheel bag (Black)	\$109.00
XP130	Frame protectors full length	\$100.00
XA150	Plastic tube fitted to offensive wing	\$71.00
XS070	Foot strap (select if required)	No Cost
XS060	Calf strap (select if required)	No Cost
XP400	Back bar protector	No Cost
XR090	Tool kit	No Cost

Spec Notes

Off Spec Requests (Please Specify)

CODE

DESCRIPTION

MSRP



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: cscanada@sunmed.com www.sunrisemedical.ca
Specifications subject to change without notice. © 2026 Sunrise Medical Canada Inc. CN-100025-CAN REV 12.