




Customer				Ship To			
Account				Name			
Account #			Date	Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City			Province
Buyer				PC	Phone		
Mark For				Email / Fax			

Please submit orders to: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)

Seating Specifications			
Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options		
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23"   Depth: 13-24"	STF Height Range	19", 20", 21"
Drive Base Specifications			
Drive Wheel Position	4 Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	400 lbs.   350 lbs. with Elevate or Tilt	Drive Motors	4 x 700W 4-pole motors
XT4 Base Size	Width 27.5"   Length 39"   Turning Radius 55.5"	Speed	6.2mph   3.1mph at full seat elevation
Min chair weight (No power options)	408lbs.	Tires	All Terrain Tires
Std chair weight (Lift, Tilt, Recline)	485lbs.	Charger	Dual Mode
Battery Size	Group 31	Ground Clearance	5"

**Client Measurements**

a) Top of head to seat	_____	in
b) Top of shoulder to seat	_____	in
c) Underarm to seat	_____	in
d) Elbow to seat	_____	in
e) Top of leg to seat	_____	in
f) Shoulder width	_____	in
g) Chest width	_____	in
i) Hip width	_____	in
j) Actual sitting depth	_____	in
m) Foot length	_____	in
n) Knee to foot	_____	in
o) Seat width	_____	in



**Notes**

✓ OPTION	DESCRIPTION	RETAIL
<b>Seating Options</b>		
<input type="checkbox"/> 600SF4	<b>REHAB</b> Seat and Back - Back Canes and adjustable flat seat pan	NC
<input type="checkbox"/> 600SF5	<b>MPS</b> Seat and Back - Captain style seating	NC
<input type="checkbox"/> 600SF6	<b>MPS</b> Back and <b>REHAB</b> Seat - Captain style backrest and adjustable flat seat pan	NC

**Base**

XT4 Base End Item Number: EIPW60 \$

All XT4 bases include four 700 watt 4 pole motors, LED lights and indicators, and ISO-7176 approved transit brackets as standard.

<input type="checkbox"/> 600MW96	Magic Mobility - XT4 No Power Seating Options	\$ 29,000.00
<input type="checkbox"/> 600MW97	Magic Mobility - XT4 Single Power Seating Options	\$ 29,000.00
<input type="checkbox"/> 600MW98	Magic Mobility - XT4 Multi Power Seating Options	\$ 29,000.00

<input type="checkbox"/> 600MW99	Magic Mobility - XT4 No Power Seating Options - HD	\$ 29,000.00
<input type="checkbox"/> 600MW100	Magic Mobility - XT4 Single Power Seating Options - HD	\$ 29,000.00
<input type="checkbox"/> 600MW101	Magic Mobility - XT4 Multi Power Seating Options - HD	\$ 29,000.00

<input type="checkbox"/> 600MW102	Magic Mobility - XT4 No Power Seating Options - MPS (Captain)	\$ 29,000.00
<input type="checkbox"/> 600MW103	Magic Mobility - XT4 Single Power Seating Options - MPS (Captain)	\$ 29,000.00
<input type="checkbox"/> 600MW104	Magic Mobility - XT4 Multi Power Seating Options - MPS (Captain)	\$ 29,000.00

<input type="checkbox"/> 600MW105	Magic Mobility - XT4 No Power Seating Options - HD MPS (Captain)	\$ 29,000.00
<input type="checkbox"/> 600MW106	Magic Mobility - XT4 Single Power Seating Options - HD MPS (Captain)	\$ 29,000.00
<input type="checkbox"/> 600MW107	Magic Mobility - XT4 Multi Power Seating Options - HD MPS (Captain)	\$ 29,000.00

**Base Accessories**

<input type="checkbox"/> 6008A91	Jack for Drive Tire	\$ 272.00
----------------------------------	---------------------	-----------

**Optional Drive Tire Options**

<input type="checkbox"/> 6008A101	Tire Sealant for 4 Pneumatic Drive Tires	\$ 103.00
-----------------------------------	--	-----------

**Shroud Color** \$

<input type="checkbox"/> 600C185	Midnight Black	NC	<input type="checkbox"/> 600C207	Digital Camo	\$ 815.00
<input type="checkbox"/> 600C188	Horizon Blue	NC	<input type="checkbox"/> 600C199	Camouflage	\$ 815.00
<input type="checkbox"/> 600C187	Sunset Red	NC	<input type="checkbox"/> 600C206	Canadian Flag	\$ 815.00
			<input type="checkbox"/> 600C209	Carbon Fiber	\$ 815.00

**Seat To Floor Height** \$

XT4/XT2 Seat to Floor Heights start at 18.75" +/- 0.25"

<input type="checkbox"/> 600HT44	Lowest Seat To Floor Height	NC
<input type="checkbox"/> 600HT45	Seat To Floor Height +1"	NC
<input type="checkbox"/> 600HT46	Seat To Floor Height +2"	NC

**Batteries** \$

<input type="checkbox"/> 600AG46	Group 31 Gel Battery (Pair) - Installed (Available with XT4)	\$ 1,450.00
----------------------------------	--	-------------

**Battery Charger** \$

<input type="checkbox"/> 600BJ38	12A Battery Charger	NC
----------------------------------	---------------------	----

### Seating: Seat Functions

✓ OPTION	DESCRIPTION	RETAIL
No Power Seat Options		\$
<input type="checkbox"/> 600ST87	No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, or future power seat functions	\$ 9,328.00
Power Tilt		\$
Utilized without Power Recline.		
<input type="checkbox"/> 600PA64	50° Power Tilt in Space - 350 lb. weight capacity	\$ 10,236.00
Power Seat Elevate		\$
Only available when Power Tilt or Power Tilt and Recline is selected.		
<input type="checkbox"/> 600PI16	12" Adjustable Seat Elevator - 350 lb. weight capacity	\$ 8,917.00
Power Recline		\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.		
<input type="checkbox"/> 600PA65	150° Power Recline Only - 300 lb. weight capacity	\$ 7,286.00
Power Tilt and Recline		\$
<input type="checkbox"/> 600PA66	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 16,165.00

### Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)	RETAIL
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.	
<input type="checkbox"/> 600SI24	16" x 18" Vinyl Seat NC
<input type="checkbox"/> 600SI25	18" x 18" Vinyl Seat NC
<input type="checkbox"/> 600SI26	18" x 20" Vinyl Seat NC
<input type="checkbox"/> 600SI27	20" x 18" Vinyl Seat NC
<input type="checkbox"/> 600SI28	20" x 20" Vinyl Seat NC
<input type="checkbox"/> 600SI29	22" x 20" Vinyl Seat \$ 617.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes	
\$	

#### REHAB Seat Pan

Standard seat pan available 13"-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 621.00
<input type="checkbox"/> 600PM15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 621.00

#### REHAB Seat Width

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

#### REHAB Seat Depth

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

#### REHAB Seating Cushion Options

Cushion size will be determined by Rehab Seat size selected above	RETAIL
<input type="checkbox"/> 600CU59	JAY® Ion® (Sizes 16" - 21" Wide) \$ 657.00
<input type="checkbox"/> 600CU60	JAY® Ion® (Size 22" Wide) \$ 698.00
<input type="checkbox"/> 600CU61	JAY® Union® (Sizes 16" - 21" Wide) \$ 764.00
<input type="checkbox"/> 600CU62	JAY® Union® (Size 22" Wide) \$ 806.00
<input type="checkbox"/> 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide) \$ 1,045.00
<input type="checkbox"/> 600CU65	JAY® Fusion® (Size 22" Wide) \$ 1,188.00

**Seating: Seat Sizes**

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Support Pads and Hardware (must pick pad and mount)		\$
---	--	----

**Lateral Thigh Support Pads**

<input type="checkbox"/> 600KF22	Lateral Thigh Support Pad Short - Right (5x5)	\$ 75.00
<input type="checkbox"/> 600KF23	Lateral Thigh Support Pad Short - Left (5x5)	\$ 75.00

**Lateral Thigh Support Mounts**

<input type="checkbox"/> 600KF26	Lateral Thigh Removable Mount - Right	\$ 325.00
<input type="checkbox"/> 600KF27	Lateral Thigh Removable Mount - Left	\$ 325.00

**Seating: Seat Back Options**

MPS Backrest Options - Comes standard with Manual Recline		\$
---	--	----

**Backrest Width**

<input type="checkbox"/> 600BN46	Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47	Wide Backrest Width - 22-23" Wide	\$ 452.00

**Backrest Height**

<input type="checkbox"/> 600B80	Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81	Tall Backrest Height - 26" Tall	\$ 744.00

**Backrest Material**

<input type="checkbox"/> 600BU106	Cloth Backrest	NC
<input type="checkbox"/> 600BU107	Vinyl Backrest	\$ 410.00

MPS Back Accessories		\$
----------------------	--	----

<input type="checkbox"/> 600PH33	Push Handle - Needed for O2 holder	\$ 211.00
<input type="checkbox"/> 6009A13	Accessory Mount - Needed for aftermarket headrest	\$ 159.00

REHAB Seating Backrest Options		\$
--------------------------------	--	----

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

**JAY J3 Back Contour**

<input type="checkbox"/> 600JP01	J3@ PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 912.00
<input type="checkbox"/> 600JP02	J3@ PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,941.00
<input type="checkbox"/> 600JP11	J3@ PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 986.00
<input type="checkbox"/> 600JP13	J3@ PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,068.00

**JAY@ J3@ Back Shell Options and Height Ranges**

<input type="checkbox"/> 600JR41	15" tall shell (Upper Thoracic Med) for 22" wide	NC
<input type="checkbox"/> 600JR14	16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15	18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR18	18" tall shell (Shoulder Height) for 22" wide	NC
<input type="checkbox"/> 600JR17	21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27	22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option		\$
------------------------------------	--	----

<input type="checkbox"/> 600OD23	REHAB Flip Forward Backrest Option	\$ 214.00
----------------------------------	------------------------------------	-----------

### Seating: Seat Back Options

✓ OPTION	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only		\$
<b>JAY® Lateral Thoracic Support Pads</b>		
<input type="checkbox"/> 600JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 75.00
<input type="checkbox"/> 600JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 75.00
<input type="checkbox"/> 600JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 75.00
<input type="checkbox"/> 600JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 75.00
<b>JAY® Lateral Thoracic Swing-Away Mount</b>		
<input type="checkbox"/> 600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 334.00
<input type="checkbox"/> 600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 334.00
JAY® Anterior Thoracic Support		\$
<input type="checkbox"/> 600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 227.00
<input type="checkbox"/> 600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 227.00
<input type="checkbox"/> 600K10	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 322.00
<input type="checkbox"/> 600K11	JAY® Anterior Thoracic Support - Large ATS8303	\$ 322.00
Pelvic Positioning Options - Select One		\$
<b>Pelvic Positioning</b>		
<input type="checkbox"/> 600A848	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/> 600A849	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 103.00
<input type="checkbox"/> 600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 159.00
<input type="checkbox"/> 600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 272.00
Headrest Options - MPS Style		\$
<input type="checkbox"/> 600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
<input type="checkbox"/> 600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
Headrest Options - REHAB Style		\$
<b>Headrest Pad</b>		
<input type="checkbox"/> 600HE39	Whitmyer® Plush Headrest Pad, 10"	\$ 377.00
<input type="checkbox"/> 600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 445.00
<b>Headrest Mount</b>		
<input type="checkbox"/> 600MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$ 247.00
<input type="checkbox"/> 600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 298.00
Armrest Style Options		\$
<input type="checkbox"/> 600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,788.00
<input type="checkbox"/> 600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,539.00
<input type="checkbox"/> 600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,167.00
Armrest Pad Options		\$
<input type="checkbox"/> 600AB81	10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE81	10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/> 600AB82	14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE82	14" Waterfall Armrest Pad - Left	NC

Legrest Options			
<input checked="" type="checkbox"/> OPTION		DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
<input type="checkbox"/> 600H210		Center Mount Flip-Up Footplate	\$ 610.00
Center Mount Fixed Legrest - Accessories			
<input type="checkbox"/> 600DI59		Fixed Center Mount Calf Pads	\$ 697.00
<input type="checkbox"/> 600LS16		Foot Strap - 1.5" wide strap	\$ 120.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest (Includes Calf Pads)			
<input type="checkbox"/> 600H215		Center Mount Power Elevating - 1 Piece footplate - Includes calf pads (11" x 11.5")	\$ 5,583.00
<input type="checkbox"/> 600H216		Center Mount Power Elevating - 2 Piece footplate - Includes calf pads (5" x 10.5")	\$ 5,417.00
Swing-Away Legrest			\$
Swing-Away Legrest			
<input type="checkbox"/> 600H211		Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 877.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
<input type="checkbox"/> 600H214		Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 6,094.00
S/A Power Legrest Positioning Options			
<input type="checkbox"/> 600HL15		Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 171.00

Accessories			\$
<input type="checkbox"/> 600A93		Waterproof Wheelchair Cover	\$ 222.00
<input type="checkbox"/> 600A94		Waterproof Joystick Cover	\$ 69.00
<input type="checkbox"/> 600A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 51.00
<input type="checkbox"/> 600A221		Deluxe Backpack	\$ 144.00
<input type="checkbox"/> 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 287.00
<input type="checkbox"/> 600A719		Articulating Oxygen Holder	\$ 1,492.00
<input type="checkbox"/> 600A612		Cane Holder - Includes push handles	\$ 957.00
<input type="checkbox"/> 6009A30		Luggage Rack for XT - Attaches to base behind the chair (10.5" x 13")	\$ 499.00
<input type="checkbox"/> 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 803.00
<input type="checkbox"/> 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 445.00
<input type="checkbox"/> 600A614		Fishing Rod Holder	\$ 445.00
<input type="checkbox"/> 6008A73		Armrest Bag	\$ 110.00
<input type="checkbox"/> 6009A10		Utility Hooks	\$ 85.00
<input type="checkbox"/> 6008A59		Transfer Bars	\$ 525.00

Be Adaptive Products			
<input type="checkbox"/> 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,242.00
<input type="checkbox"/> 6008A99		Be Adaptive Trigger Mechanism	\$ 1,242.00
<input type="checkbox"/> 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 2,242.00

Electronics			
Joystick			\$
<input type="checkbox"/> 600J121		Standard LED Joystick (w/o Screen Display)	NC
<input type="checkbox"/> 600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
<input type="checkbox"/> 600J118		High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$ 757.00

Electronics											
✓ OPTION		DESCRIPTION						RETAIL			
<b>Joystick Mount</b>											
\$											
<input type="checkbox"/>	600JM34	Fixed on Armrest - Depth adjustable 8.25"						NC			
<input type="checkbox"/>	600JM45	Swing-Away Joystick Mount - Includes latch to lock in the forward position						\$ 849.00			
<input type="checkbox"/>	600JM46	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)						\$ 2,050.00			
<b>Joystick Mount Location</b>											
\$											
<input type="checkbox"/>	600HM16	Right Hand Armrest						NC			
<input type="checkbox"/>	600HM17	Left Hand Armrest						NC			
<b>Display Language Selection R-net</b>											
\$											
<input type="checkbox"/>	600L1	English	NC	<input type="checkbox"/>	600L3	French	NC	<input type="checkbox"/>	600L7	Spanish	NC
<b>Joystick Handle Type</b>											
\$											
Comes standard with PG Cone-Shaped Joystick Handle											
<input type="checkbox"/>	600JH13	Ball Handle						\$ 110.00			
<input type="checkbox"/>	600JH26	Bodypoint - Grooved Mushroom Joystick Handle						\$ 240.00			
<input type="checkbox"/>	600JH24	Bodypoint - Rubber Dome Joystick Handle						\$ 256.00			
<input type="checkbox"/>	600JH25	Bodypoint - U-Shaped Joystick Handle (4 inch)						\$ 256.00			
<b>Expandable Power Module</b>											
\$											
Must select both items											
<input type="checkbox"/>	600EX05	Expandable 120A Controller - R-net						\$ 1,099.00			
<input type="checkbox"/>	600EX06	Harness for 120A Expandable Controller						\$ 768.00			
<b>Power Seat Function</b>											
\$											
<input type="checkbox"/>	600EO37	Thru-Drive Control for 1 Actuator <small>Must be selected if ordering 1 actuator function</small>						\$ 1,443.00			
<input type="checkbox"/>	600EO38	Thru-Drive Control for 2+ Actuators <small>Must be selected if ordering 2 or more actuator functions</small>						\$ 2,903.00			
<b>Power Seat Control Option</b>											
\$											
<input type="checkbox"/>	600QG11	Ctrl+5 Switch Box with Free Assignable Functions						\$ 1,017.00			
<input type="checkbox"/>	600QG12	Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)						NC			
<b>Ctrl+5 Switch Type</b>											
\$											
<input type="checkbox"/>	600EY22	Push Button						NC			
<input type="checkbox"/>	600EY23	Toggle Switches						\$ 356.00			
<b>Power Seat Control Option - Ctrl+5 Mounting location</b>											
\$											
<input type="checkbox"/>	600EY19	Right Armrest	NC	<input type="checkbox"/>	600EY20	Left Armrest	NC				
<b>Attendant Controls</b>											
\$											
<input type="checkbox"/>	600A05	Dual Control Module for Attendant - Includes accessory mount and push handle						\$ 2,460.00			
<b>Attendant Control Mounting Location</b>											
\$											
<input type="checkbox"/>	600DN5	Right Back Cane	NC	<input type="checkbox"/>	600DN6	Left Back Cane	NC				
<b>R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared &amp; Bluetooth</b>											
\$											
<input type="checkbox"/>	600ON14	R-net Omni-2 Only Input Control Module						\$ 3,611.00			
<input type="checkbox"/>	600ON15	R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices						NC			
<b>R-net Omni-2 -Mounting Location</b>											
\$											
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC	<input type="checkbox"/>	600DS24	Left Hand Armrest	NC				
<b>Head Control Packages</b>											
\$											
<input type="checkbox"/>	600EU25	Switch-It DualPro® Proportional Head Array - Facial Pads						\$ 7,057.00			
<input type="checkbox"/>	600EU26	Switch-It DualPro® Proportional Head Array - Spot Pads						\$ 7,057.00			