



## Clinical Seating designed for kids!

The JAY Zip® cushion and backrest are specially designed to meet the unique sizes and needs of kids.

The dual layered, contoured foam base of the **ZIP CUSHION** provides skin protection and positioning in a comfortable, effective, and convenient package.

The **ZIP BACKREST** is an ultra lightweight backrest featuring antimicrobial cover technologies and the adjustability of J3 Back RH hardware.



**zip**

**SUPERIOR CLINICAL SEATING**

## Verstile, Quick-Release Hardware

- Ultra lightweight and durable mounting hardware
- Fits on a wide range of tubing, from 3/4" to 1 1/8"
- Up to 2.75" of depth adjustment
- +/- 10° of angle adjustment

## Convenient and Antimicrobial

### Outer Cover

- Silver impregnated X-Static® stretch fabric naturally inhibits bacterial growth
- Extra layer of microclimatic spacer fabric for comfort
- Machine washable and quick to dry
- Two outer covers included (cushion only)
- Safe non-toxic materials

### Inner Cover

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base

## Comfortable

- Dual layered, contoured foam base
- Reticulated foam in the outer cover promotes air movement
- Soft, stretchable outer cover

### JAY TECHNOLOGY IN THE ZIP CUSHION AND BACK



## JAY ZIP SPECS

### Cushion

HCPCS Code  
E2607

Weight Capacity  
165 lbs

Weight  
1.2 lbs (12x12)

Width Range  
8" to 16"

Depth Range  
8" to 18"

Height  
2.75"

### Back

HCPCS Code  
E2613

Weight Capacity  
165 lbs

Weight  
2.5 lbs (12x12)

Available Widths  
8", 10", 12", 14"  
Adjustable +1"

Available Heights  
6", 8", 10", 12"

Depth Adjustability  
2.75" (anterior)

Back Angle  
+10° to -10°



©03.2019 Sunrise Medical (US) LLC  
MK-129528 Rev.A

800.333.4000    www.SunriseMedical.com

