



**ORDER FORM**

Customer Dealer Information		Shipping Address	Identical To Billing Address
Dealer		Client	
Client		Address	
Address			
E-Mail:		E-Mail:	
Phone:		Phone:	
Order Information		User Information	
Purchase Order No.:		Current Wheelchair:	
Special Requests (INTERNAL USE ONLY):		User Disability:	
		User Weight:	kg      User Height

SIGNATURE: \_\_\_\_\_

DATE: \_\_\_\_\_

**Serial No:**

**RGK Octane Sub4**

SBAXZ010

TIGA Sub4

\$9,056.00

**Max User Weight 110kg (242lbs)**

**Please Note: All RGK chairs are hand built and have a manufacturing measurement tolerance of +/- 5mm**

**Product Style**

XA070

Open Frame

No Cost

XA080

Closed Frame (Non certified Sub4 option)

No Cost

**CAD Approval**

XW009

Quote only (No CAD For Approval)

XW010

Order CAD Drawing For Approval  
*(2 Customer Amendments ONLY)*

\$300.00  
FLAT FEE

XW011

Order (No CAD Required)

XW012

Quote & CAD For Approval  
*(2 Customer Amendments ONLY)*

\$550.00  
FLAT FEE

**Key Measurements**

XY010

Seat Width (Inside of sideguards)

mm (Min - 250mm / Max - 440mm)

XY020

Seat Length

mm (Min - 250mm / Max - 440mm)

XY030

Extra Seat Length (Provides additional leg support and hand grip)

mm (Max - 440mm Total Seat)

No Cost

XY050

Front Frame Height

mm (Min - 380mm / Max - 520mm)

XY060

Rear Frame Height

mm (Min - 360mm / Max - 500mm)

XY070

Frame Length

mm (Min = XY020 + XY030 +140mm / Max - 625mm)

XY072

Castor Position



mm Optimise Measurements

XA020

Ergonomic Seat

mm Specify Length - (Standard = 120mm)


\$121.00



## Front Frame

XY080	<b>Seat Taper Design</b> (When combined with XY091, dimension is inbetween the centre of the main bends) mm	<b>No Cost</b>
XY091	<b>Parallel Taper Design</b> (When combined with XY080, dimension must be smaller by 20mm) mm	<b>No Cost</b>

## Footrest Measurements

XY100	<b>Round Design</b> (Castor wheel size 4", max height achievable 110mm)		<b>Std</b>
XY105	<b>Square Design</b> (Castor wheel size 4", max height achievable 130mm)		<b>No Cost</b>
XY115	<b>Front Footrest Height</b> mm		
XY120	<b>Rear Footrest Height</b> mm    STD = 20mm lower than front height		
XY125	<b>Footrest Length</b> mm    Maximum Sub4 length = 140mm (certified weight) (160mm and over requires heavy duty upgrade XB150)		

## Footrest Style


XB015	<b>Tubular Aluminium Footrest Sub4 - Height Adjustable 15.88 x 1.6mm</b>	<b>Std</b>
XB010	<b>Tubular Aluminium Footrest Non Sub4 - Height Adjustable 16 x 3mm</b>	<b>No Cost</b>
XB200	<b>High Level Footrest</b>	<b>\$201.00</b>

mm    Indicate front footrest height



## Footrest Cover

XB045	<b>Lightweight Carbon Cover</b> (With weight saving cut-out) <b>No Cut-out</b> (Non certified Sub4 Option) <b>Footplate Width Sub4 certified</b> (STD if no option selected) <b>Footplate Width Standard</b> (Optimised plate width) (Non certified Sub4 Option)	<b>Std</b>
-------	---	------------



## Castor Wheels


XC135	Aluminium Centre 4" x 1.25" Silver Hub / Black Tyre (Std) Black Hub / Black Tyre Purple Hub / Black Tyre Green Hub / Black Tyre Red Hub / Black Tyre Blue Hub / Black Tyre	Std
SBAXC155	Forged Carbon Centre 4" x 1.4" 	No Cost
XC140	Aluminium Centre 5" x 1.4" (Silver Hub/Black Tyre) 	\$232.00
XC145	Aluminium Centre 6" x 1.4" (Silver Hub/Black Tyre) 	\$310.00

## Castor Forks

XC190	Performance Castor Fork (Std) <small>Note: Colour follows accent</small>	Std
XC192	High Trail Castor Fork <small>Note: Greater distance between contact patch and steering axis reduces flutter at higher speeds. Colour follows accent</small>	No Cost

## Castor Pin Width

XY330	Castor Pin Width <small>Optimised Measurements (Minimum width for clear castor rotation and maximum wheelbase length)</small>
-------	--



## Wheel Lock Options

SBAXD036	Compact Scissor with Lightweight Fixing	No Cost
XD020	Underseat Scissor 	No Cost
XD050	Push to Lock 	No Cost
SBAXD085	Eagle Wheel Locks 	\$220.00
SBAXD086	Eagle Wheel Locks with L/W Fixing  <small>(Only with fixed centre of gravity)</small>	\$240.00
	Black (Std)	
	Purple	
	Red	
	Blue	
XD010	Omit Wheel Locks	

## Sideguards

XE055	Screw on Lightweight Carbon Level with Tyre Individual Above Tyre (Follow Tyre - Max 25mm) mm		No Cost
XE065	Screw on Lightweight Carbon with Fender (Fender centre of tyre) 10mm Above Tyre (Std) Individual Above Tyre (Min 10mm) mm		No Cost

## Backrest

XG010	Fixed Backrest (Non Folding)		Std
SBAXG020	Folding Backrest with 2 Position Snap Lock & Angle Adjust		\$716.00
XY230	Backrest Height (Max 400mm) mm		
XY360	Backrest Angle -5 -3 0 +3 +5 Individual (Max -10/+20 Deg) Deg		
XG100	Ergonomic Backrest Vertical Height mm Rearward Position mm		\$146.00
XG110	Inward Taper (Inside Of Backrest Posts - Specify Width Or Per Side) mm		\$146.00
XG120	Outward Flare (Inside Of Backrest Posts - Specify Width Or Per Side) mm		\$146.00

## Camber Bar

XY390	Camber Angle 1 2 3		
SBAXH100	Single Post Fixed Centre of Gravity (25-125mm) mm		Std
SBAXH130	Square Back Adjustable Centre of Gravity (25-125mm) mm		\$866.00
SBAXH140	Round Back Adjustable Centre of Gravity (60-125mm) mm		\$866.00

## Push Handles

XG070	Screw in Push Handles		\$225.00
XG080	Flip Down Push Handles		\$441.00

## Anti-Tip

XJ100	Single Sport Anti-Tip		\$504.00
-------	-----------------------	---	----------

## Frame Finish

SBAXN030	Brushed Frame	No Cost
SBAXN050	Polished Frame	No Cost
XN010	Powder Coat Frame 	No Cost

Paint Code

## Accent Colour

Accent colour includes, forks & castor caps

Black  
Red  
Purple  
Silver  
Blue

## Seating Options

XP045	Sub4 Seat Sling Strap Adjustable	Std
XR050	Carbon Fibre Seat Plate  mm (Std XY020 - 20mm)	\$480.00

## Sub4 Backrest Upholstery Options

XP035	Airtech Sub4 backrest strap adjustable	Std
-------	--	-----

## Non Sub4 Backrest Upholstery Options

XP050	Airtech backrest strap adjustable 	No Cost
XP070	Backrest upgrade to leather corners (Black)  (Only with airtech upholstery XP050)	\$121.00
XP110	Backrest with zip pocket  (Only with airtech upholstery XP050)	\$127.00
XP170	Individual backrest embroidery (Please supply JPEG)  (Only with airtech upholstery XP050)	\$240.00

## Backrest Upholstery Stitching Colour

Silver	Red
Black	Orange
Blue	Purple
Individual	

## RGK Logo Colour

White (Std)	Red
Black	Orange
Blue	Purple
Individual	

### Rear Wheel Size

XY420	24" Wheels	XY430	25" Wheels
XY460	Rear Wheel Gap -	15	20 25

### Rear Wheel Style

XK075	Spinergy BLXL 18 Blade Spoke (N/A in Silver Hub)	Std
SBAXK150	Spinergy CLX 18 Carbon Blade Spoke (N/A in Silver Hub)	\$1552.00
SBAXK040	Supply Without Wheels	Wheel Type: Tyre Type
XK226	Titanium 100mm quick release axle pins	

### Spinergy Spoke Colour

	Black	Red	Purple	
	White	Blue	Orange	
	Yellow	Pink	Green	
XX180	Individual Spoke Design	Inside	Outside	\$213.00

### Spinergy Hub / Rim Colour

XK190	Black Hub (Std)	Std
XK200	Silver Hub	\$127.00
XK201	Red Hubs & Rims	\$409.00
XK202	Purple Hubs & Rims	\$409.00
XK203	Blue Hubs & Rims	\$409.00
XK204	Gold Hubs & Rims	\$539.00

### Handrims Style

XL010	Aluminium hard anodised 6 tab	Std
XL160	Ellipse 6 Tab	\$631.00
XL020	Tetra Grip 6 Tab	\$400.00
XL100	Para Grip 6 Tab	\$400.00
XL040	Titanium 6 Tab	\$446.00
XL035	Natural Fit LT	\$743.00
XL170	Carbolife Curve L (BLK)	\$1393.00
XL150	Carbolife Gekko - Black Silver	\$248.00

### Tyre Style

XM010	Kenda - Grey Blue Red (Grey standard if no option selected)	Std
XM060	Schwalbe One (Lightest Option)	\$133.00
XM020	Schwalbe Marathon Plus	\$98.00
XM040	SHOX Solid (26" N/A with Fusion Wheels)	\$160.00
XM130	Kenda Black	\$28.00
XM050	Kenda K885 mountain bike tyre (24",25" only)	\$40.00

### Inner Tube

XY500	Schraeder (Car Valve) - Standard if no option selected	 XY500	No Cost
XY490	Presta (Sports Valve)	 XY490	No Cost

### Wheel Build Up Options/Accessories

XY480	Narrow Set Handrim		
XY485	Wide Set Handrim		
XL090	Remove Excess Mounting Tabs	\$146.00	
Accessories:			
XK145	Black Gloss Spokeguards	Fitted Supply Non - Fitted	\$400.00

## Cushion

XP080	Cushion Foam Soft (CMS040) Medium (CMS050) 50mm 75mm 100mm	\$121.00
-------	---	----------

## Cushion Cover

XP085	Nylon	Std
XP090	Airlite	\$98.00
XP100	Leather	\$254.00

## Extra Frame Options

XR060	Welded calf strap mounting loops <small>(Specify Height from Floor)</small> <small>Optimise Measurement</small>	\$154.00
-------	--	----------



## Accessories

XR080	Frontwheel Assembly (Only with Tubular End Angle Adjustable Footrest) - Select Colour: Black (Standard) Blue Silver Purple Red	\$1,053.00
XP120	Frame protectors Neoprene zip fasten (Black)	\$84.00
XP420	Frame protectors Nylon velcro fasten (Black)	\$81.00
XR010	Seat pouch velcro attach (Black)	\$59.00
XR200	Rucksack with RGK logo (Black)	\$151.00
XR090	Tool Kit	No Cost
XS065	Sub4 Calf Strap	No Cost

Spec Notes

Off Spec Requests (Please Specify)

CODE

DESCRIPTION

MSRP



1000 Creditstone Road, Unit 2, Concord ON L4K 4P8

Customer Service: (800) 263-3390 Fax: (800) 561-5834 E-mail: [cscanada@sunmed.com](mailto:cscanada@sunmed.com) [www.sunrisemedical.ca](http://www.sunrisemedical.ca)  
Specifications subject to change without notice. © 2026 Sunrise Medical Canada Inc. CN-100038-CAN REV (12). Page 9 of 9