



PERSONALIZE YOUR JAY[®] MADE TO ORDER SEATING

Vinyl is a water resistant, easy to clean surface. It is recommended for planar seats and backs. Vinyl can be used on non-contact surfaces to add color to your seating system.


Stretch upholstery is recommended for contoured seats and backs, anti-thrust seats and similar items. **Darlexx[®]** and **Dartex[®]** allow four-way stretch that is anti-microbial, anti-fungal, waterproof and breathable. Dartex is two-sided allowing you to choose between a cloth surface or a smooth surface that is easy to wipe clean. **Neoprene**, sometimes referred to as “wet-suit material”, allows four-way stretch. Neoprene is a soft nylon material with rubber backing. It is waterproof but not breathable.

Cordora is a heavy-duty, textured nylon fabric that dries quickly if wetted.


Nylon is a finer weave and smoother texture than Cordora but offers the same excellent durability.

3DX[™] Microclimatic protects the skin by dissipating heat and moisture.

Stretch




Cordora





Nylon



3DX[™]



Texture Vinyl



See page 2 for other Vinyl options

Vinyl

